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3.2.9.	Temperatures of the cycle.....	45
3.2.10.	Timer of the cycle	47
3.2.11.	Hospital data and serial number	49
3.2.12.	Inverter Ammeter Setup.....	49
3.2.13.	Hertz Settings.....	51
3.2.14.	Inverter Scale.....	52
3.2.15.	Cycle Counter	53
3.2.16.	Temperature Sensor Setup	53
3.2.17.	Bypass Inspection door.....	54
3.2.18.	Manual Controls.....	55
3.3.	Switching on and switching off	55
3.4.	Automatic cycle.....	56
3.5.	Manual machine controls	57
3.6.	Service messages.....	58
3.7.	Alarms.....	59
3.8.	Filling of sterilization vessel	63
3.9.	Carry out a pseudo-cycle in 'Manual'	64
3.10.	Downloading material from the cell in MANUAL mode.....	66
3.11.	Verifying the effectiveness	67
3.12.	Integrated Waste Collection system.....	68
4.	Maintenance.....	71
4.1.	Informative notes.....	72
4.2.	General norms.....	72
4.3.	Replacing consumable material.....	74
	List and duration of consumable parts.....	74
4.3.1.	List and duration of consumable parts.....	74
4.3.2.	Changing printer paper.....	79
4.3.3.	Vertical blade changing.....	80
4.3.4.	Dismantle blade support.....	81
4.3.5.	Changing blades.....	82
4.3.6.	Replacing short temperature sensor.....	82
4.3.7.	Replacing the short sensor cover.....	84
4.3.8.	Replacing the long temperature sensor.....	85
4.3.9.	Replacing the long temperature sensor cover.....	85
4.3.10.	Changing the lid gasket.....	86
4.3.11.	Replacing the carbon and absolute filter.....	87
4.4.	At the end of the cycle.....	88
4.5.	End of the work session.....	88
4.6.	Maintenance for long inactivity.....	89
5.	Problem solving.....	91
5.1.	Rotor blocked because of overload.....	92

5.2.	The printer does not print and the cycle does non finish.....	93
5.3.	The treated material is too hot when discharged.....	93
5.4.	Strange noise of metal banging inside the cell.....	94
5.5.	The machine heats up less than usual	94
5.6.	Steam escapes from the cell.....	95
5.7.	Table of the main problems, their causes and solutions.	96

1. GENERAL INFORMATION

1.1. Identification of main parts


Fig. 1 Layout


Fig. 2 Layout - Side A


Fig. 3 Layout - Side B


Fig. 4 Vessel group


Fig. 5- Demister group


Fig. 4 - Filter Group

Fig. 5 - IWC- Integrated waste collector

1.2. Identification plate

For any communication to the manufacturer or the distributor, it is necessary to refer all the information on the identification plate.

The identification plate, including the machine serial number, is located on the top of the machine, in a visible position.

1.3. Use of the machine

The purpose of the NEWSTER 5 is to treat potentially infectious solid healthcare waste in order to obtain:

- Sterilization
- Physical modification
- Dehydration
- Volume reduction
- Weight reduction

The machine and its components must be used only for the treatment of solid healthcare waste, containing liquids, potentially contaminated by pathogenous microorganisms, deriving from hospitals, laboratories, medical research centers, dialysis centres, and infectious disease hospital wards as specified below [according to DPR n.254- 15 of July 2003- Attachment I (art.2, comma1, lettera a)]

European code	Kind of waste	Classification
CER 180103 or 180202 Wastes coming from research activity and bacteriological diagnostic	Sanitary towels, napkins and nappy Cotton flock for papetest and colposcopy Ocular sticks not sterile TNT ophthalmic sticks Cannulaes and drain pipes Catheters (vesical, venous, arterial, pleural drains, pipe-unions, probes) Extrabody circulation circuits Expired blood bags Defluxors Contaminate phleboclysis Dialysis filters, exhausts air filters (without chemical risk) Surgery gloves Disposables : vials, pipettes, tubes, protective dresses, mask, glasses, lenghts of material, sheets, foot wear, white overalls Materials for medication (bandages, tampons, wax plasters, lunghette) Bags (transfusionals, nutritionals, urine)	Dangerous with infective risk

European code	Kind of waste	Classification
	Infusion sets Rectal and gastric probes Probes (nose, broncoaspiration, oxygen therapy) Auricular speculums Vaginal speculums Plaster casts or bandages Teeth and small body parts not recognizable Bedstead for experimental animals Empty containers Wastes coming from dentists Wastes from restaurant (infective only) Petri plates, and all contaminated disposables	
Cutting wastes CER 180103 or 180202	Needles, syringes, knives, venflon, razors and bistoury disposables	Dangerous with infective risk
Cutting wastes not used CER 180101 or 180201	Needles, syringes, knives and bistoury	Special
Anatomic wastes CER 180103 or 18202	Tissues, organs and anatomic parts not recognizable Animals for experiments	Dangerous with infective risk
Empty containers CER	Empty containers from drugs, veterinary drugs, disinfectants, food, drinks and infusion solutions	Special/assimilable to civil wastes

The machine is not suitable for treating:

- Gas bottles or similar
- Chemical products
- Inflamable material
- Explosive material
- Radioactive material
- Stones
- Bulky or hard pieces larger than 100 gr.
- Any material similar to the above listed

1.4. Environmental conditions and operating limits

The machine can be installed in an ordinary room, even of small size, provided with aeration, an electrical outlet of adequate power and a connection to the water grid and to the sewage.

To ensure proper operation the sterilizer unit must be installed in an area sheltered from the weather, with an ambient operating temperature between 5 ° C and 45 ° C and a relative humidity not exceeding 95%. The machine can work 24 hours per day.

1.5. Noise Level

Noise measurements were carried out according to UNI EN ISO 3746 and UNI EN ISO 11202.

The noise level was found to be of: 70 db (A)

1.6. Technical specification

Thermal sterilizing treatment:	Proteinic demolition in humid environment after dehydration, with temperatures as high as 150° C at the end of the cycles
Nominal waste disposal potential:	15 kg/ hour of waste with 10% humidity 240 liters/hour
Final volume of treated waste:	20-25% of initial volume
Final weight of treated waste:	70-75% of initial weight
Overall control system:	Programmable logic controller (PLC)
Discharging temperature	90 °C
Water consumption:	About 50 liters/day with the water recirculation system
Temperature measurement:	By variable-resistance sensors
Cycle recording:	Recording of temperature and time
Final cooling of waste:	At approximately 100° C through evaporation of

	H ₂ O
Dust removal:	In humid environment with separation on cartridge filter
Cell lid closure:	Manual with mechanical block and double safety device
Safety devices:	Mechanical blocking of sterilization vessel lid. Direct interruption of main motor power supply in case of failed lid block. Low voltage commands with automatic power cut-off if panels are opened.
Transmission of movement to rotor:	By pulleys and V-drive belts
Rotating blades and fixed counterblades:	In special steel
Water discharge:	Through tap to sewers with 40 mm diameter pipe
Electrical power requirements:	Maximum overall requirements: 20 kW.
Dimensions and weight:	
Sterilizer L x P x H	160 x 80 x 130 cm Kg 580
Filter storage L x P x H	80 x 35 x 150 cm Kg 60
Electrical board L x P x H	75 x 45 x 110 cm Kg 100
Weight:	740 kg

Emissions:

Air emissions discharged from the machine conform with the environmental protection requirements, as witnessed by the certification located at the headquarters of the manufacturer.

Before being discharged, the air is treated by active carbon filters and finally by an absolute filter. The absolute filter has an DOP >99,99% efficiency.

The water that is discharged in the sewer conforms to the requirements, as witnessed by the certification located at the headquarters of the manufacturer.

Noise level is cut down by polyurethane shields.

The machine does not produce any waste other than the waste treated.

1.7. Disuse of the machine

When the machine is not in use, all the machines components must be disposed of conforming to the rules and regulations of the country where the machine is installed. The recyclable parts must be separated and disposed of at the proper disposing facility. Electrical components must be separated from the rest of the material.

1.8. Warranty

For warranty conditions see the attached form that must be filled and sent to Newster System Srl.

2. SAFETY NOTES

2.1. Preliminary observation

Access to the rooms assigned for the use of the machine must be permitted only to personnel authorized by the healthcare management authority.

You are reminded that in the rooms assigned for use of the machine it is forbidden to drink, smoke or eat.

For all operations required, the use of protective gloves is compulsory.

The various risk situations that might affect the machine have been analyzed, verifying the precautions and technical solutions required to eliminate or reduce them in the most appropriate way by applying adequate measures of protection.

2.2. Transport and movement of machine

With an empty weight of 580 kg and overall size and shape envisaged, various types of lifting equipment can be easily used. The machine is factory-assembled and packaged in suitable dimensions for transport. The machine can be positioned on any normal floor capable of supporting the relative loads, and four adjustable feet are provided to compensate for unevenness of the surface.

2.3. Proper use of the sterilizer "Newster® NW05"

The device and its components may only be used for the treatment of solid potentially infectious hospital waste as listed in paragraph 1.3.

Do not use the equipment to treat gas containers, chemical, flammable, explosive or radioactive materials. Hard and compact masses (eg. compact metal bodies with mass greater than 100 grams, stones, etc..) can cause damages to the rotating and fixed blades.

The apparatus is not suitable for the treatment of liquids or waste having less than 60% of dry matter. In these cases, the waste may be treated after being mixed with material from a previous treatment.

The apparatus can not be loaded only with materials with a melting point of less than 150°C (as PVC). In this case the waste can be treated by mixing with another material which has a melting point higher than 150°C. In any case, the average weight of materials with a melting point less than this value should not exceed 60% of the total weight.

The equipment is designed to operate at atmospheric pressure (in order to enable the cell to be in depression). Do not cover, for any reason, the exhaust outcoming vapours, otherwise the system may go slightly under pressure.

2.4. Qualified Personell

Operating on this equipment is allowed only to trained personnel, namely:

- operational staff who has received proper training to work with the equipment and who knows the contents of this manual and of the annexed documentation,
- service and maintenance personnel shall be trained to repair automation equipments and authorized to connect and disconnect from electrical power, open, ground and work on electrical circuits, components and systems according to what is established by safety regulations.

Do not leave the equipment available to persons who are not properly trained.

The non-compliance of the information contained in the manual or affixed on the equipment itself can cause serious damage to people and property.

After opening the panels of the machines and / or of the electrical board panels, certain parts of the machine become available, they can have dangerous

voltages, which could result in death or serious injury to persons and property if appropriate precautions wouldn't be taken.

2.5. Protection against accidents

The device is designed with the following protection against accidents. They should not be removed for any reason:

- electrical equipment follow the standards EN 60204-1
- the auxiliary circuits are made in low voltage (24 V),
- sterilization of the cell operates under a pressure lower than atmospheric .
- the main electrical disconnecting-switch and electrical board door are interlocked,
- an electrical interlock prevents the starting of the rotor if the lid is not closed properly,
- an electrical interlock prevents the opening of the lid if the rotor is in motion,
- to close the lid, the operator must first unlock a safety piston,
- the system provides a self-check to verify the closure of the discharging door,
- in case of lack of tension, the lid of the cell and the discharging door remain tightly closed,
- in case of failure the machine stops in the best conditions of hygiene and safety allowed by the situation.

2.6. Labels and safety symbols

Labels bearing symbols and/or safety messages are applied on the sterilizing machine. In this section we show their exact location and the kind of hazard.

Attention. Read the user manual carefully before using the machine or doing maintenance.

Figure 8

Danger. Moving parts. Danger of crushing arms during operation of the main motor. Never put your hands under the sterilizing machine while the main engine is working.

Figure 9

Figure 10

Danger. Risk of crushing upper limbs during the inspection and cleaning of the discharge mouth. Never put your hands inside the inspection hole of the discharge during the closing of the mouth. Close the inspection hole by acting on the outlet with both hands. To clean the drain use the tool provided with the machine (see figure 11).

Figure 11

Danger. Hot surface. Danger of burns of the upper limbs during the operation of the machine. Do not hold the lid by pulling on the metal parts, open and close the lid acting only on the plastic handles.


Figure 12

Attention. Make sure that all the labels are always readable. Wipe them with a cloth with soap and water. Do not use solvents. Replace the damaged plates receiving them from the manufacturer. If a plate is located on a piece that must be replaced ensure that the new one is of the same type and is located in the same position.

2.7. Residual risks

The careful analysis carried out by the manufacturer and stored in the file, made it possible to eliminate most of the risks associated with the conditions of use of the steriliser. The manufacturer recommends to strictly follow the instructions in this manual and safety regulations, including the use of protection devices required.

The residual risks can be:

Electrical hazards during maintenance under tension

Following the DLgs 17/2010 only trained personnel is allowed to operate under tension and in accordance with the double insulation to ground.

Risk of improper use in the presence of explosive and flammable substances

Follow the instructions in paragraph 2.3.

Risk of operating in a wet environment

Comply with the limits given in section 1.4.

Risk of crushing of the upper limbs during the cleaning of the discharge

The machine has a label indicating the security risk. It is possible to avoid this risk by using the special tool as shown in section 2.6.

Risk due to moving parts

The machine has a label indicating the security risk. It is possible to avoid this risk by acting as specified in paragraph 2.6.

Risk of burns of the upper limbs

The machine has a label indicating the security risk. It is possible to avoid this risk by using the special plastic handles as indicated in paragraph 2.6.

Risk of misinterpretation of the pictograms

The symbols are installed in accordance with the rules for graphic symbols to be used.

The user must immediately replace the safety signs that may become unreadable due to wear or damage.

Attention. It is absolutely forbidden to remove the safety labels on the machine. The manufacturer declines all responsibility in the event of breach of this prohibition.

The machine described in this manual applies to the following guidelines:

2006/42/EC, 2004/42/EC and 2006/95/EC.

2.8. Authorized personnel and responsibility

To prevent injury, this apparatus must be managed, by the person responsible, one operator at a time.

When the equipment is turned on, the support staff and maintenance must remember not to put their hands under the edges of the lid of the cell and must perform operations only after obtaining the consent of the responsible conductor at that time.

2.9. Hygiene norms

Always remember that the handling of hospital waste involves the potential of incidental contact with pathogenic microorganisms. So: at work wear the protection clothes (to be replaced at least weekly and washed with disinfectant) and use disposable gloves when handling bags containing waste. If the wastes are originated from infectious wards, increase caution during these operations, by also wearing a mask, a pair of glasses and respiratory protection.

Keep the working space well ventilated (at least 10 air changes per hour).

Do not smoke, keep or consume food in the workplace.

Operators and the technical assistance should not touch their face, especially the mouth, nose and eyes before washing hands with disinfectant soap.

Always wash your hands with disinfectant soap after work and before eating, drinking, smoking (in a separate room).

Even though before the end of each cycle the equipment sterilizes itself, for caution before performing maintenance on parts which can come into contact with the waste or vapors, sprinkle abundantly with disinfectant (disinfectants used: based on chlorine, iodine, hydrogen peroxide and the like), close the lid and turn on the heating elements, wait 15 minutes for the disinfectant to act, turn on the vacuum cleaner filters and after 5 minutes turn everything off and then begin to maintenance work. At the end of maintenance work the protective coat needs to be changed.

In case you have an accidental spill of liquid on the floor, proceed to the collection and cleaning with disinfectant.

At the end of the work, wash the floor with water and disinfectant.

2.10. Protection from electric shock

The equipment must be connected to a good grip on the ground before being turned on.

The operation and maintenance of electrical components and automation equipment may be carried out only by trained personnel.

The power cables to the board are always in tension. For eventual maintenance to the input clamps of the electric board the operator must remember to interrupt the power supply to the electric board

Never throw water on live electrical components.

Always turn off the machine at the end of the working day and in any case before maintenance you must make sure that the electrical circuit is completely de-energised by

turning the disconnecter, located on the left of the front panel, horizontally to 'O' (Or 'OFF').

2.11. Emergency stop

In the case of urgent need to stop the device, press the button **EMERGENCY STOP** obtaining the immediate arrest of all the components, for the fall of the auxiliary circuits.


Figure 13

Figure 14

The buttons "EMERGENCY STOP" are two and they are: in the lower right corner on the console of the machine (Figure 13) and on the control panel as shown in (Figure 14).

2.12. Electric welding on the equipment

Electric welding on the structure or any part of the equipment can cause serious damage to electronic components. Request for service and disconnect the PC from the electrical panel.

2.13. Operations on components sensitive to static electricity

Do not allow unqualified personnel to do maintenance to electronic components. The basic rule to follow is that electronic components should be touched with the hands only if you need to work to be performed on them. Do not touch the components contacts and their conductors.

Touch components only if your body is electrically discharged continuously to the ground by means of a wrist strip conductor or special shoes or shoes with antistatic conductive strips to the ground.

Before touching an electronic module, be sure to not be electrostatically charged. The easiest way is to touch something grounded (such as a water pipe) before touching the module.

2.14. Preventing fires

Keep a dry powder fire extinguisher on hand. Water extinguishers are not suitable.

Never spray water on electrical components.

To avoid risk of combustions, before leaving the workplace unattended, ensure that the treated material is not hot over 100 ° C.

If the material is discharged at temperatures above 130 ° C, it must cool with water and it must absolutely not be close in the collection bag until its temperature is not less than 80 ° C. A material treated with too high temperature could cause an a fire.

If the material being processed has a strong smell of roasted, it means that it has been subjected to excessive temperatures (over 180 ° C) due to malfunction. The most likely causes are the lack of water during the cooling phase, or anomalies to the system of temperature measurement. Control the transport system of the water and the chain of temperature regulation.

2.15. "Logical" equipment

The normal operation of the equipment is in "Automatic" mode via a logic and a sequence of phases and interlocks that follow one after another. This logic is factory set and is contained in the EPROM of the CPU of the PLC. The times can

be changed only by specialized programmers (if authorized and provided the password) and are held in RAM memory by the lithium battery. If you remove the battery to restart it will load the default time data stored in the EPROM from the manufacturer.

Changes to the electronic logic involves a risk of malfunction and may only be performed by experienced programmers in automation.

2.16. Locked rotor stress

No particular caution must be taken during the normal operation of the machine because its interior is accessible only when the system has passed through a state of sterility. However, in case of locking without possibility of restarting before the completion of the sterilization cycle, it is necessary to use some caution.

Wear a lab coat, heavy-duty gloves (eg. Leather) and a mask with filter (follow the hygiene rules in this manual) when you open the cover always be protected by gloves, glasses and mask and empty the machine by pushing the material through the exhaust port in a few bags that can later be reworked into a new cycle.

If the failure occurred immediately after the start of the cycle, and if the material is not yet sufficiently shredded to be extracted through the discharge port, it is unfortunately necessary to proceed manually from above. Wear a lab coat, heavy-duty gloves (eg. Leather) and a mask with filter (follow the hygiene rules in this manual).

2.17. Operating precautions

All control devices are clearly visible, and the operating cycle has been planned to simplify and limit the activity of the operator to the turning on of the machine,

the loading of the waste into the sterilization vessel and the starting of the treatment cycle.

To avoid all risks, operation in automatic mode (cycle selector set to automatic) for the entire cycle provides a suitable guarantee. However, it is also possible to select manual mode for the cycle to check operation without load or to solve problems that arise during automatic operation.

The risk of access to the operating zone of the machine (vessel with rotor and blades) has been eliminated by:

- A system that locks out start of cycle if the lid is not closed (with circuit microswitch);
- A system that locks out start of cycle if the discharge port hatch is not closed (with circuit microswitch);
- A system that locks the lid in place during the machine operation;
- A system that blocks the opening of the discharge hatch during operation of machine.

The potential risk of breakage of the rotating blades or parts of the rotor support system (rotating flange and blade holder) are analyzed under normal operating conditions.

The components of the rotor are sized according to the stress envisaged for the treatment by pulverization of materials currently used in hospitals (in particular needles, scalpels, sanitary boxes, metallic probes, glass, as well as plastics, cellulose products, organic wastes, etc). The rotor is therefore capable of resisting the impacts and stress foreseen for its operation.

However, to avoid damage, every time the blades are removed the fixing bolts must be replaced. The replacement bolts must be of the same type.

The machine is equipped with a button for normal stopping and a suitable emergency stop button with a manual reset system to return the machine to operational status.

The following possible malfunctions of the machine during operation were considered:

- In case of a power failure, the lid and the discharge valves for materials remain locked in place to prevent their opening or the emission of materials from the vessel when the machine cycle starts again.
- In case of partial malfunctions after the start of the cycle (lack of pressure in water circulation system, partial electrical malfunction, failure to register temperature, etc) the machine is stopped and the problem is warned of on a specific display, so that it can be rectified before the automatic cycle is restarted.
- The control and monitoring system is easy to use, with easily identifiable buttons, simple instructions and alarm messages shown on display or with indicator lights, with ample explanations in the manual.

The electrical system of the machine is compliant with the relative standard, harmonized in the EU zone. The machine must be used only when connected to a ground cable and system, as envisaged in the installation phase (inside the door of the electrical panel). The power supply line must be protected with a delayed-action circuit breaker set to trip above 100 A. Power supply cables for the machine must have a cross-section of 10 mm² if cable length is shorter than 20 meters, 16 mm² if cable length is between 20 and 50 meters, and 20 mm² if cable length is over 50 meters. Electrical installations must be carried out by a qualified expert in electrical systems.

At the end of daily operating sessions the electrical system of the machine must be turned off with the main switch located on the electrical panel (position "0"), and this operation must normally be carried out by the user of the machine.

To avoid the risk of spontaneous ignition, at the end of daily operating sessions it is essential to check that any treated materials remaining in the vessel are at a temperature of less than 100°C.

When handling healthcare wastes, which create the hazard of contact with pathogenic micro-organisms, the operator must wear protective overalls and gloves to handle the bags containing wastes. Overalls must be made of a fabric resistant to acids and alkalis, and must be changed at least once weekly and

sent for hot washing in chlorine bleach. Gloves (compliant with EN 374) must be of the disposable type. If wastes from isolation wards are to be treated, when the vessel lid must be opened or there are wastes or residues inside the vessel during internal vessel cleaning operations, all handling and movement operations must be carried out wearing a face mask (compliant with EN 140) to protect the respiratory tract, in addition to overalls and gloves.

The room in which the machine is installed must be ventilated with an extractor that performs at least 10 exchanges per hour of the total air volume. The room must also be equipped with a washbasin, fire extinguisher and running water. The operator must be warned to refrain from touching the face (mouth, nose and eyes) unless the hands have been washed first. At the end of every work session it is always compulsory to wash the hands with a disinfectant soap with povidone-iodine (Betadine or similar product). Provide warning signs and forbid smoking and the keeping of food in the room.

2.18. Maintenance precaution

For maintenance operations, and in particular for maintenance operations inside the vessel, the main switch must always be turned off.

To deactivate the electrical system the main switch on the electrical panel must be moved to the "0" position.

Maintenance operations on the electrical system must be carried out by specialized personnel.

Electrical panel is not interlocked with the opening lever of the panel, that means the power supply terminals are still live if the electrical panel is opened.

All maintenance operations must be carried out using the personal safety equipment envisaged for this type of work. At the end of maintenance operations, protective overalls must be changed.

A dry-powder extinguisher with a capacity of 6 kg must be kept in the vicinity of the electrical panel of the machine or at some other accessible point in the room where the machine is installed, for use in case of possible emergencies.

2.19. List of required precautions and equipment

Personal safety equipment (to be worn in its entirety during maintenance operations on parts in direct contact with wastes):

- protective overalls resistant to acids and alkalis (to be hot-washed in chlorine bleach once weekly);
- disposable gloves (compliant with EN 374);
- face mask (for cleaning vessel);
- protective glasses (for handling NaClO) and eyewash.


2.20. Other precautions:

- washbasin with running water and disinfectant soap with povidone-iodine (Betadine or similar products); availability of a canister of water;
- dry-powder extinguisher capacity 6 kg;
- warning signs (including prohibition of smoking and keeping of food in the room);
- information to operator on hygiene precautions to be taken (for mouth, nose, eyes);
- ventilation of the room (10 exchanges of air volume per hour).

3. INSTRUCTIONS FOR THE USER

3.1. Description of process

The NEWSTER®NW5 converter-sterilizer is designed to treat Healthcare Risk Waste (HCRW) at atmospheric pressure and high temperature in a wet environment, as illustrated in the following diagram:

- First stage: loading and start of heating phase
- Second stage: evaporation of liquids
- Third stage: superheating and sterilization
- Fourth stage: cooling
- Fifth stag : unloading

In a closed sterilization vessel a powerful rotor fitted with blades disintegrates, agitates and heats the wastes by impact, friction and resistances.


The temperature of the mass of wastes is measured in real time and with great accuracy by special sensors, while they are being constantly and vigorously agitated.

When the temperature reaches the predetermined level of 150°C, the mass of wastes is automatically sprayed with water so that the temperature decreases.

The treated wastes are cooled down to 95°C. At this point the cycle has been completed, and the product, by now sterile, is automatically unloaded.

To disperse the heat produced by the system, part of the water is continuously replaced by fresh water from the mains supply. Excess water and incondensable gases are discharged into the sewers, with values within legally established limits.

The process lasts for about 25 minutes, depending on the rate of organic matter in the waste, and automatically follows these stages:

1. Wastes are loaded into the sterilization vessel, the lid is closed and the treatment process is started by pressing a button. The engine starts in the first speed, filters start and resistances are turned on.
2. At 60°C the general water electrical valve opens and water starts to flow in the cooling columns. The rotor turns slowly at first, starting to pulverize the materials, and at the same time the temperature starts to rise.
3. The rotor turns faster (change from the first speed to the second), the temperature starts to rise rapidly and the materials are finely pulverized.
4. When 96–100°C are reached, the temperature remains stable until the water present in the wastes has completely evaporated.

5. After the water has evaporated the temperature starts to rise rapidly again, reaching 150°C. After the peak the rotor turns slowly (first speed) and resistances turn off.
6. The mass of wastes is sufficiently moistened by a water spray to cool down to 95°C.
7. The sterilization cycle has now been completed. The vessel is opened, and the product is extracted and collected in the stainless integrated waste collector.

3.2. Comands


Fig 3.1

3.2.1. The operator panel

The normal operation of the equipment will take place in “automatic” mode by means of a logical linking of phases and blocks that follow each other in an automatic sequence.

The logic is pre-loaded at the factory and is contained in the programmable logic controller (PLC) on an EPROM card. This card contains pre-loaded default operating parameters.

Only an authorised user with the key software is able to enter into a dialogue with the system and to introduce changes to the parameters.

WARNING! Modifications to the electronic logic carry the risk of disfunction and thus making the guarantee void. In any case, such modifications must only be made by programmers who are experts in automation.

The operating panel permits the control and variation of the timers and temperature values which control the normal course of the automatic cycle. It can also send alarm messages and information on the state of the machine.

The operating panel is touch screen (figure a.)

Figure a.

You can check the status of the cycle at any time by touching the icon.


Figure b.

3.2.2. Main menu

By pressing the icon  at any time the main menu will appear.


By pressing the icon  at any time the alarm messages menu will appear as described in the proper chapter.

Figure c.

From the main menu you can enter 4 submenus.

Icon  to display the inner cell temperature (par. 3.2.3)

Icon  to display the system submenu (language, date and time, input/output etc) (par. 3.2.4)

Icon  to display the layout (par.3.2.1)

Icon  to display the cycle parameters set up in the submenu (par.3.2.7)

Icon  to display the manual controls buttons (par.3.2.18)

3.2.3. Display inner cell temperature

Use this page every time you want to know the temperature inside the cell.

Temperatures below 40°C will not be shown.


Figure d.

Pressing the function button you can go to 2 submenus

Icon  to display the cycle parameters set up in the submenu (par.3.2.7)


Icon  to switch manually the motor in first speed.


3.2.4. Sub menu system

Figure e.

Instructions for the user

Icon  to set date and time (par. 3.2.6).

Icon  to change the language (par. 3.2.5).

Icon  to change the panel contrast (par. 3.2.7).

Icon  to see the PLC INPUT and OUTPUT state.


Figure f.

Other function buttons are used only by expert technicians. Refer to the user manual of the PLC KTP700 Siemens to find out about these functions.

If you have questions regarding the use of the products described in the manual and are unable to find an answer here, please contact your nearest Newster representative

3.2.5. Language selection

Language selection

Figure g.


Select the language by pressing the function button beside the language chosen or the icon. PLC changes the language immediately.

3.2.6. Set up time and date


Figure h.

Touching the blank field you want to fill out this window will appear:

Touch again on the blank field an alphanumeric keyboard appears. Use it to write all the data.

Press the icon  when finished, and to start the clock.

Press the icon  to go back to settings submenu.

3.2.7. Display Contrast Setup

Figure i.

Press the icon of the desired display contrast value to set the contrast.

Press the icon  to calibrate the touch screen

3.2.8. Set up cycle parameters sub menu



Figure I.


Icon  to set up temperatures of the cycle submenu (par. 3.2.9).

Icon  to set up the timer of the cycle submenu (par. 3.2.10)


Icon  to set up hospital data and serial number (par.3.2.11)

Icon  to set up the current value for speed change (par.3.2.12).

Icon  to set up temperature sensor (par.3.2.16)

Icon  to set up cycles number and motor hours (par.3.2.15).

Icon  to set up inspection hatch microswitch (par.3.2.17)

Icon  to display the inverter hertz set up (par.3.2.13)

Icon  to display the inverter scale set up (par.3.2.14)

3.2.9. Temperatures of the cycle



Figure m.


Icon  to display the cycle parameters set up in the submenu (par.3.2.7)

In these submenus you can set up the temperatures of the cycle, furthermore these temperatures are printed to describe the course of the cycle. Not all values are possible, inside PLC are configured some limits and every temperature must respect these limits.

It is highly recommended to contact the producer before changing these values to avoid any malfunction during the cycle.

Severe personal injury or substantial property damage can result if proper values are not entered in this data base.

Touching the blank field you want to fill out this window will appear:

Touching the blank field again an alphanumeric keyboard appears. Use it to write all the data.

In this data base it is possible to set the following values.

- Opening general water: general water valve opens at this temperature and water starts to run inside heat exchanger column.
- Pump restart: pump with deodorizer (where it is present) starts again at this temperature.
- Sterilization temperature: it is the maximum temperature reached during the cycle
- Closing water in cell: at this temperature the discharging process starts during the cycle
- Max temp external cell wall: It is the higher temperature that can be reached on the external wall of the cell. At this temperature the cycle will stop and an alarm will appear
- External temperature of the cell: It is the temperature on the external wall of the cell
- Start water cooling vapours: at this temperature cooling water will start to flow
- Stop water cooling vapours: at this temperature cooling water will stop to flow

3.2.10. Timer of the cycle

Figure n.



Figure o.

Icon  or  to move between the 2 set up timers of the cycle submenus.


Icon  to display the cycle parameters set up in the submenu (par.3.2.7).

In these submenus you can set up all the timers of the cycle. Not all values are possible, inside PLC some limits are configured and every timer must respect these limits.

It is highly recommended to contact the producer before changing these values to avoid any malfunction during the cycle.

Severe personal injury or substantial property damage can result if proper values are not entered in this data base.

Touching the blank field you want to fill out this window will appear:

Touch the blank field again an alphanumeric keyboard appears. Use it to write all the data.

In this data base it is possible to set up the following values:

- Main motor in first speed: after this time the motor changes speed.
- Time ON pump first start: at the beginning of the cycle, a pump with deodorizer (when present) works for the time indicated.
- Time ON pump second start: During the cycle a pump with deodorizer (when present) starts again and works for the time indicated.
- Temperature delay OK: once that temperature is measured on the PLC, the temperature changes after the time indicated, in order to prevent that a noise signal has any effect on the cycle.
- Discharge door delayed opening: when the discharging temperature is reached PLC waits the time indicated before opening the hatch, to prevent vapour leaks.
- Product discharge time: the discharging phase will last the time indicated.

3.2.11. Hospital data and serial number




Figure p.

Icon  to display the cycle parameters set up in the submenu (par.3.2.7).

In these submenus you can set up all the hospital data and the serial number of the sterilizer. This data is printed on the reports.

The hospital data must stay on 2 lines of 20 characters each; meanwhile the serial number can have only 3 numbers.

Touching the blank field you want to fill out this window will appear:

Touch the blank field again and an alphanumeric keyboard appears. Use it to write all the data.

3.2.12. Inverter Ammeter Setup



Figure q1.

In these submenus you can set up the current's value for motor speed change. Touching the blank field you want to fill out this window will appear:

Touch the blank field again and an alphanumeric keyboard appears. Use it to write all the data.

3.2.13. Hertz Settings



Fig. q2

Touching the blank field you want to fill out this window will appear:

Touch again on the blank field an alphanumeric keyboard appears. Use it to write all the data.

Press the icon  to go back to cycle parameters settings submenu.

3.2.14. Inverter Scale



Fig. q3

Touching the blank field you want to fill out this window will appear:

Touch again on the blank field an alphanumeric keyboard appears. Use it to write all the data.

Press the icon  to go back to cycle parameters settings submenu.

3.2.15. Cycle Counter

Figure r.

In these submenus you can see the cycle counter and the working time of the motor.

3.2.16. Temperature Sensor Setup



Figure s.

In these submenus you can set up the sensor parameters. Not all values are possible, inside PLC some limits are configured and every value must respect these limits.

It is highly recommended to contact the producer before changing these values to avoid any malfunction during the cycle. Severe personal injury or substantial property damage can result if proper values are not entered in this data base.

Touching the blank field you want to fill out this window will appear:


Touch the blank field again an alphanumeric keyboard appears. Use it to write all the data.

3.2.17. Bypass Inspection door

Figure u.

Touching the slide  icon this window will appear:

Touch the blank field again and an alphanumeric keyboard appears. Use it to write username and password. Now it is possible to activate or deactivate the safe alarm for inspection door and for the recirculating system.

We warn that is possible to bypass this alarm only for maintenance by authorized operator.

3.2.18. Manual Controls


Figure v.

In these submenus you can activate in manual some functions illustrated in par.

3.5

3.3. Switching on and switching off

Put the mode selector on "Manual" and turn the "Main Switch" knob to position 1. Check for presence of power supply (power light indicator on).

At the beginning the CPU loads the internal program in 30 seconds, at the end the display shows this message:

EMERGENCY PRESSED OR MACHINE NOT STARTED

Verify that emergency stop buttons are not pressed, otherwise put them in the normal position rotating the red button, push the button "START MACHINE" (the button green light must switch on), now the machine is started.

The display will show again the same message above, only after that the "START MACHINE" green light is on press the button ACK on the operator panel, the message will disappear.

If after this the machine will not start check the solving problem section.

In case of emergency machine can be switched off by pushing one of the red emergency stop buttons.

Warning! Use this option only in emergency.

To switch off the machine press the "STOP MACHINE" button and turn the main switch on "0" or OFF.

Warning! During normal operation the "BY-PASS SAFETY DEVICES" lamp must be off, otherwise call authorized personnel in charge of the maintenance.

3.4. Automatic cycle

Put a bag under the discharging hatch in the steel collector and close it. With the cell full and the lid closed rotate the key selector to reach "Automatic" mode, wait 10 seconds until the machine closes the hatch (if the user forgot it open) and

check the safety devices. Press the button . The automatic cycle will start.

If the machine does not start, check the display on the PLC panel for error messages by pressing F2. Check the display and follow the instruction in paragraphs 3.6 and 0.

Cycle can begin only if the display shows

NO MESSAGE

The signal for the end of the operating cycle is given by a red light located above the main control panel, which gives both a light and sound signal. To halt this signal, turn the mode key to reach "Manual" position. When the cycle is finished

and the bagging machine has completely stopped (where this accessory is present, before taking out the bag it is important to shake the filter of the bagging machine so all the material attached to the filter can drop inside the bag). Now the machine is ready for another cycle.

WARNING! Always wear anti-dust mask to close and change the sterilized material bags.

3.5. Manual machine controls

To use the 'Manual' functioning mode, all you have to do is rotate the selection mode key on "MANUAL" position.

By pressing buttons on the keyboard and on the manual functions submenu (par. 3.2.18) the following functions can be carried out:

Start the automatic cycle


LED on - start motor at low speed, LED off - STOP motor


LED on - filters fan in operation, LED off - ventilator stopped

LED on - resistance operating, LED off - resistance off

LED on - general water valve open, LED off - valve closed

LED on - water valve in cell open, LED off - valve closed

LED on - door in opening phase, LED off - door closed

LED on - door in closing phase, LED off - door closed

LED on - lid unlocked, LED off - lid locked

Instructions for the user

Increases the speed motor

Decreases the speed motor

LED on - motor turns in clockwise, LED off - motor turns in counterclockwise

LED on - pump is working, LED off - pump stops

Note that when the LED is green it means that the function is on in manual mode, while when the LED is blue it this means that the function is on in automatic mode.

The normal operation of the machine will be in 'Automatic' mode.

The operations carried out in 'Manual' mode are exclusively for

- checking the functioning of the system or of other components when empty or when loaded,
- solving problems which occur during the execution of the automatic cycle,
- carrying out emptying operations and cleaning the sterilization cell.

3.6. Service messages

To access the level of alarm signalling, press the icon . In this way you enter the alarms menu and by scrolling through the messages using the arrows you will find that any alarms present are shown on the display.

If the display shows this message.

HANDLE OPEN

It means that the handle is open, in this case the automatic cycle cannot start

and in manual mode the motor cannot rotate. Push the key  (verify the key

has the led on), close the handle and press the button again, this time the led must be off, force the handle a little. If the message is still showing, contact authorized technical assistance.

If the display shows this message.

LID OPEN

In this situation the automatic cycle cannot start and in manual mode the motor cannot rotate because the cover is not closed. Even in this case push the key

and verify that the led is on, take the handle of the cover and close it. The message should disappear.

If the message is still on contact authorized technical assistance.

If the display shows this message:

DISCHARGING HATCH NOT BLOCKED

The discharging hatch is not closed correctly, the automatic cycle cannot start.

Put the machine in manual mode, open the discharging hatch pushing

clean the surface of the hatch well and close again with the key

The message should disappear, otherwise contact the technical assistance.

When pushing ESC the display shows this message.

NO MESSAGES

Now the cycle can start normally.

3.7. Alarms

In this paragraph the main alarm messages and the procedure for resetting the machine for normal operation after an alarm message are examined.

To access the level of alarm signalling, press the icon . In this way you enter the alarms menu and by scrolling through the pages using the arrows on the screen you will find that any alarms present are shown on the display.

If during the normal operation of an automatic cycle the PLC detects a reason for an alarm, the cycle will immediately be interrupted and the display will begin to flash indicating the alarm, a red light will go on in the alarm column and you will also hear an acoustic signal. To reset the correct functioning of the machine, follow one of the next procedures, according to the nature of the alarm.

LIST OF ALARMS

Immediately after turning on the general power switch and after having pressed EMERGENCY STOP buttons, the screen displays the following message.

EMERGENCY PRESSED OR MACHINE NOT STARTED

Verify that the two “**EMERGENCY STOP**” buttons are not pressed, otherwise replace them to the normal position by triggering them off with counterclockwise turn.

Press “**START MACHINE**”, the button should light up, now press the ACK icon and the message will disappear.

If you press the **START MACHINE** button and the light does not turn on verify that the “**EMERGENCY STOP**” button is not pressed and check that the protection panels of the machine are tightly screwed in. Press **START MACHINE**, if the machine still does not start then call technical assistance.

If the display shows

THERMAL PROTECTION FUMES SUCTION

turn off the power, positioning the main switch in the “0” (off) position, open the electric board and find the breaker marked MT2 and move its lever to the “1” (on) position. Then close the switchboard and turn the main switch in the “1” (on) position. Finally, follow the instructions at the bottom of this list of alarms.

If the display shows

MAXIMUM CURRENT MAIN MOTOR PROTECTION

turn off the power, positioning the main switch in the "0" (off) position, open the electric board and find the breaker marked MT5 and move its lever to the "1" (on) position. Then close the electric board and turn the main switch to the "1" (on) position. Finally, follow the instructions at the bottom of this list of alarms.

If the display shows

THERMAL PROTECTION FILTER SUCTION MOTOR

turn off the power, positioning the main switch in the "0" (off) position, open the electric board and find the breaker marked MT1 and move its lever to the "1" (on) position. Then close the electric board and turn the main switch to the "1" (on) position. Finally, follow the instructions at the bottom of this list of alarms.

If the display shows

LOW WATER PRESSURE

check that both the water tap is open and that there are no obstructions along the hose that carries the water to the machine. Then follow the instructions at the bottom of this list of alarms.

If the display shows

OVERLOAD INVERTER

it means that a fault stops the drive. You must reset the drive and find a solution to the problem.

Before you contact the distributor or the factory because of unusual operation, prepare some data. Write down all the texts on the display, the fault code, the fault ID, the source info, the Active Faults list and the Fault History.

To reset the drive follow the instructions at the bottom of this list of alarms.

TEMPERATURE SENSOR BROKEN

it means that in 22 minutes the temperature in the cell do not rise up to 40 °C. If you are at the first cycle of the day and that the machine was cold it is possible that the sensor is not damaged and it is only necessary to restart the cycle.

If the problem remains contact the assistance service.

If the display shows

If display shows

CELL TEMPERATURE TOO HIGH DO NOT OPEN THE LID

it means that the cell's internal temperature is excessively high, at a level which could damage its external covering. In this case you need to wait until the temperature shown on the indicator falls to 90°C before switching on the machine again (follow the instructions at the bottom of this list of alarms).

WARNING. For no reason open the lid and pour in water in manually.

Warning! This alarm could be due to a malfunction in the temperature probe. Before using the machine again in the normal way we advise that you consult the appropriate paragraph in the chapter regarding maintenance and, if necessary, contact the assistance service.

If display shows:

INSPECTION DOOR OPEN

check the inspection door and if it is necessary, close them. Then follow the instructions at the bottom of this list of alarms.

Once the cause of the alarm has been eliminated, the display continues to flash until the alarm has been acknowledged. To acknowledge the alarm, press icon

. If the LED is lit because an alarm is indicated, pressing this key will switch it off, and cause the panel cursor to stop flashing.

However, if other alarms are also present, the display will continue to flash to indicate them, in which case repeat the sequence just described until all have been removed.

If after this procedure the display shows the alarm again, please contact authorized technical assistance.

3.8. Filling of sterilization vessel

Put the machine in manual mode, open the hatch by pressing  and clean the metallic surface around the discharge door with the special device supplied

with the machine, then close the hatch again by pressing the button , now

press the button  (button led lights up). With this operation the lid handle

will be released. Open the lid and then press the button  again (this time the led must be off) so the lid is blocked, now fill the vessel with the potentially infectious waste.

When filling the vessel, the lighter bags must be placed above the blades to facilitate the pulverization process, meanwhile the heavier must be placed last. If possible fill the cell with material from different wards.

The vessel must be completely filled but without exerting pressure on waste, do not compress waste.

Please remember that the machine is for solid waste and the cycle cannot finish with humidity higher than 40%.

After the vessel has been filled with wastes, press  (the led must be on) and close the lid and handle, now press  again (the LED must be off) and force the handle forward. Now the machine is ready for a cycle.

DO NOT FILL THE VESSEL WITH THE DISCHARGING HATCH OPEN.

WARNING! THE LID SURFACE CAN BE VERY WARM. GRIP THE LID ONLY USING THE PLASTIC HANDLES.

DO NOT FILL THE VESSEL WITH THE LID UNBLOCKED. THE LED ON THE KEY  MUST BE OFF.

WARNING! ALWAYS WEAR AN ANTI-DUST MASK TO CLOSE AND CHANGE THE STERILIZED MATERIAL BAGS.

3.9. Carry out a pseudo-cycle in 'Manual'

Although, as previously mentioned in the present manual, the normal functioning of the machine is in 'Automatic' mode, it is also possible to carry out a pseudo-cycle in 'Manual' mode to check the machine or finish a failed cycle.

To do this, first of all load the material to be treated. Make sure that the loading door is closed, and secure the lid tightly.

Put the machine in manual mode.

Start the rotor by pressing the  key, and at the same time start up the fan of filters by pressing , the resistance  and the pump .


After 60 seconds, press the  key once again. Now the green LED will be off and the pump will stop working.




When the temperature reaches 50 °C, press .

Observe the ammeter. Once the consumption of electric current has fallen below 100 A, wait 2 minutes while checking that the ammeter does not exceed this threshold. If it does not, increase the speed motor by pressing the  key until the inverter gains 5 Hz.

If the ammeter does not exceed 29 A repeat this operation after 2 minutes.

Observe the ammeter. If the ammeter does exceed 34 A decrease the speed motor by pressing the  key until the inverter loses 5 Hz.

Leave the rotor to rotate until the temperature of 150°C is reached. At this point press the  key until in the operator panel you can read 20 Hz and the  key in order to switch off the resistance.


After having reduced the motor speed, press the  button and cool the mass until the temperature falls to about 95 °C. Press the  button again in order to stop the cooling and wait a minute in order to allow the water to evaporate. After stopping the "water in cell", please check if "general water" is activated, otherwise activate it by pressing .

At this point it is possible to unload the treated material by pressing the  button.

The unloading of the material will take about 2 minutes.


Once the unloading has finished, stop the rotor by pressing the  button.

Allow 2 minutes to pass, then close the unloading hatch by pressing the  key.

Now press the  button and operate the pump for 30 seconds.

This operation causes disinfectant to be sprayed into the cell.

Wait 1 minute to let the spray be completely sucked into the absorption column,

then press  and  to switch off the ventilator and the electrical valve that pours water in cell.

From this point on it is possible to open the lid and start loading operations for a new cycle.

3.10. Downloading material from the cell in MANUAL mode.

Sometimes the automatic cycle may be interrupted after reaching a temperature of 150° C, in this case it is not necessary to begin the cycle again as the material is already sterile, it is sufficient to empty the vessel using the following procedure.

Make sure the internal temperature of the vessel is inferior to 95° C, that a bin with an empty bag is placed under the discharge hatch and the lid for discharge

inspection is in position. Start the motor at slow speed by pressing ,

set the speed at 20 Hz pressing the  button and after 10 seconds open the

discharge hatch , wait 2 minutes then press the  button, the motor should turn off . Now the cell is empty you can start new cycle.

Hz.

Never discharge material with the motor on speed higher than 20

NEVER OPEN THE LID TO INSPECT MATERIAL WHEN THE MOTOR IS ON, FRAGMENTS OF STERILE WASTE MAY BE VIOLENTLY EJECTED

3.11. Verifying the effectiveness

The sterilization cycle of the machine can be verified by using biological indicators.

Inside the vessel (Fig. 5) there is a compartment for the biological indicator, before filling the vessel unscrew the cap (A in Fig. 5) with a hexagonal 17 mm key and insert the testing vial. Be sure to tighten the cap properly to avoid the possible breakage of the vial.

Always use leather gloves as inside the vessel remains hot even after cycle.

Refill the vessel and perform an automatic cycle. At the end of the automatic cycle open the lid and using leather gloves and an hexagonal key unscrew the cap again and extract the vial with caution.

The vial is available for sterilization verification.

Fig. 6

WARNING! USE ONLY GLASS BIOINDICATORS WITH DIAMETER NOT MORE THAN 11 MM.

WARNING! THE NEWSTER 5 STERILIZER DOES NOT WORK UNDER PRESSURE, THEREFORE THE EFFECTIVENESS CANNOT BE VERIFIED BY USING BIOLOGICAL INDICATORS THAT ARE ACTIVATED ONLY BY PRESSURE ABOVE ATMOSPHERIC PRESSURE.

3.12. Integrated Waste Collection system

The bagging system is integrated in NEWSTER®NW05. The sterilized material is automatically discharged into a steel box connected by an airtight seal to the discharge system.

At the end of the cycle, after the discharging phase is completed, the steel box is taken away using the support with wheels.

In case of trouble with the discharging system an alarm appears on the control panel (see paragraph 3.6).

It is necessary to replace the sack for every cycle.

The automatic bagging system

The discharging sequence:

4. MAINTENANCE

4.1. Informative notes

The complete maintenance of Newster 5 machines cannot be covered in these paragraphs. A thorough maintenance can only be carried out by qualified personnel who have attended a training course either by the producers or authorized distributors.

4.2. General norms

Maintenance must be carried out by specialized technicians who have been trained in the following specific areas:

- Mechanical maintenance
- Electrical maintenance

The professionalism and competence of the technician, lies with the person responsible for safety.

Before starting maintenance intervention the person in charge of safety must:

- Vacate the work area of those not involved.
- Make sure that the necessary instruments are at hand and in good condition.
- Verify that there is sufficient light, otherwise provide portable 24 volt lights.
- Make sure the authorized personnel are equipped with the safety regulation clothes (protective gloves, goggles etc).
- Verify that the maintenance technician has carefully read the instruction manual and understands the function of the machine.

Before starting the technician must disconnect all supplies (electrical and pneumatic) and ensure the security block on the machine is in function.

Warning! If it is absolutely necessary to work with the machine during the maintenance, the service engineer must stay at a safe distance and the safety blocking switch must be close at hand.

After the intervention and before restarting the service engineer must check the security system and the complete performance.

Work on the electrical motors must be executed by ONLY by the service engineers who have been trained by either the producer or the distributor of the machine.

Warning! After every service the engineer is responsible for the complete check and safety of the machine and the protective operative system

Maximum dependability of the machine at a minimal cost is the result of a maintenance program and a careful inspection during the life of the machine

Warning! Before starting any maintenance or cleaning service the machine's energy sources must be shut down, meaning the disconnecting switch and the air compressors loading valve.

Warning! Before starting any maintenance or cleaning service the engineer must be wearing the protective clothing, goggles and gloves.

Warning! Before starting any type of service on the machine affix a sign on the machine indicating:

**MACHINE MAINTENANCE
DO NOT ADD MATERIAL**

Warning! The mechanical instruments inside the cell are very sharp, therefore always use protective gloves.

Before the re-starting operation of the machine, check the entire system in accordance to the start up program. After every service carry out a few trial cycles thus verifying the correct function of the machine.

If these safety precautions are not followed severe harm may be caused to personnel.

Before any work:

- Check the correct pictogram is displayed on the machine.
- Check that the water tap of the machine is open.
- Check the water discharge tube of the machine.

- Check if any tools were left inside the machine.

Any intervention to the sterilizer by person or persons not qualified or not authorized by producer or distributor caused the guarantee to be void.

Use original spare parts or parts advised by the producer or the distributor.

4.3. Replacing consumable material

This paragraph will examine the material that the operator needs to replace. For maintenance interventions that are not treated in the manual please contact your after sales service.

Any intervention made on the sterilizer by unqualified or unauthorised personnel will make any warranty made at the sale of the machine void.

List and duration of consumable parts

4.3.1. List and duration of consumable parts


In the following tables you will find the stimated duration of consumable parts.


Due to a stronger use of the machine it is possible a shorter duration of the parts. Please check, at the end of each work session, the consumption of the consumable parts.

GRNW 010-03-03				
CODE	DESCRIPTION	N° CYCLES	N° HOURS	MAX
MC00207	Radial Shaft Seal, AS50x70x10 VITON	600	400	
FC00222	Screw TCEI Zinc 10x20	400	267	
MC00572	Flange Upper Radial Shaft Seal	3000	2000	
MC00592	Support Blades	1500	1000	
MC00593	Support Vertical Blades (Bag Braker)	3000	2000	
MC00596	Blade Cast Iron	200	133	


MC00600	Plate Blade Support Holder	1500	1000	
MF00603	Vertical Blades	1000	667	
MC00223	OR Seal Ø126.6	3000	2000	
FX00224	Screw, TE 16x55	1000	667	
FX00268	Screw, TCEI M10x45	1000	667	
FC00188	Galvanized Nut, M10	400	267	
MC00604	Anti release plate	1000	667	
MC00211	OR Seal Ø158.11	3000	2000	
FX00283	Screw, TCEI M6x16	1000	667	
FX00284	Washer, Ø6	1000	667	

GRNW 005-04-01 UPPER CELL GROUP				
CODE	DESCRIPTION	N° CYCLES	N° HOURS	MAX
AC00098	Spray nozze, water	4000	2667	
AC00099	Spray nozze, hypochlorite	500	333	
MC00276	Seal For Lead L=1500/mt	600	400	
MC00348	Seal V-ring, V-65A	600	400	


GRNW 010-09-03				
CODE	DESCRIPTION	N° CYCLES	N° HOURS	MAX
MC00369	Absolut Filter for Newster NW 10	3000	2000	1 year
MC00368	Filter, Active Carbon	1500	1000	6 months


SANW 005-00-00

CODE	DESCRIPTION	N° CYCLES	N° HOURS	MAX
MC00157	Net, Filled Bodies Block	2000	1333	
MC00156	Filled Bodies	4000	2667	


SANW 005-03-00				
CODE	DESCRIPTION	N° CYCLES	N° HOURS	MAX
MF00667	Circular Filter For Lid	1500	1000	


GRNW 005-00-00				
CODE	DESCRIPTION	N° CYCLES	N° HOURS	MAX

Maintenance

MC00456	Sonda PT100	2000	1333	
MC00601	Fixed Blade Tempered	1500	1000	
MC00739	Cover Short Sensor, 35	250	167	
MC00751	Cover Long Sensor, 35	250	167	
MF00742	Short Sensor	2000	1333	
MC00752	Protection Measuring Hole Cell,35	250	167	

GRNW 005-06-02				
CODE	DESCRIPTION	N° CYCLES	N° HOURS	MAX
MC00616	Plate Blocking Gasket Door	1000	667	
MC00617	Seal Discharging Hatch	1500	1000	

4.3.2. Changing printer paper

Use only thermal printer paper. Width 57,mm (+/-0,5mm) and the diameter should be less than 50mm. Change paper as follows (see diagram below):

- Open the door of the printer and press on the support of the tilting mechanism of printer where it indicates PUSH;
- Insert the tip of mechanism of the press and place the roll of paper, respecting the direction of rotation as shown in figure;
- Press the FEED button to let a few centimeters out of the printer and push again on the support mechanism of printer where it indicates PUSH, then return to original position;
- Rip the paper and close the door by inserting the card slot in the door itself.

- Also take out MC00593.
- Clean MC00593 and all contact surfaces.
- Check all screws. Change if it's necessary.
- Install new vertical blades. Use only original parts.
- Place the parts in their place, tighten the screws.

Warning. Parts installed inside the cell have sharp edges. Always use protective gloves when working.

4.3.4. Dismantle blade support.

Blade support (MC00592) handles strong dynamic events. Specially for this, it is made with a special hard and resistant steel. It should be changed every 1000 cycles.


Fig. 2

- Dismount the vertical blades, as is described in chapter 4.3.3.
- Unscrew with hexagonal key on extension bolts 221.
- Raise the blade support from cell. Use protecting gloves.
- Clean all contact surfaces well before the installation of the new blade support.
- Mount all parts in their places. Tighten the screw with the correct torque.

Use only original or authorized spare parts. Using non-original or non-authorized parts may cause malfunction of the machine, or accidents.

Warning. Parts installed inside the cell have sharp edges. Always use protective gloves when working.

4.3.5. Changing blades.

Fragmentation and heating of the material is made by blades. For that reason, to keep maximum efficiency of the machine, is important to pay extra attention about the condition of the blades. Change them every 500 cycles.

- Dismount the blade support from the machine, and place it in stable vise.
- Use an allen key 8 and a size 17 wrench to unscrew all 12 bolts . Bolts are locked by nuts placed on the bottom side of blade support.
- Clean contact surfaces of the blade support and blades.
- Change used bolts with new ones.
- Place new blades and tighten the screw reliably.

222 Blades screws

Fig. 3

Warning. Parts installed inside the cell have sharp edges. Always use protecting gloves when You work.

4.3.6. Replacing short temperature sensor.

The short temperature sensor (called after short sensor) is responsible for the temperature monitoring signal which is transmitted to the PLC and printer, to control the cycle.

The machine is designed in way to protect short sensor. But it can happen, that replacing will be needed. To do so , follow the next procedure.

- Switch off the machine and take off the panels.
- Disconnect the two wires connected to the short sensor. (Fig. 4).
- With a size 7 and 17 wrench, release the locking (Fig. 5), it is not necessary to unscrew it completely, it is enough to loosen it .
- With a size 7 wrench unscrew the short sensor (Fig. 6).
- Place the new sensor, make sure that the head of the sensor touches the wall.
- Lock nut 10 with a size 17 while holding the short sensor with a size 7.
- Connect the wires to the sensor. The Way of connection is not important.
- Close the machine's panel. Make sure that the protection circuit is active.

Short sensor is one of the parts responsible for the correct work of the machine. Use only a sensor supplied by the producer of the machine, or authorized sellers.

Figura 4

Figura 5

Figura 6

4.3.7. Replacing the short sensor cover.

Because of the reason explained in chapter 4.3.6 it is very important to protect the short sensor, so please pay special attention about condition of the short sensor cover. Replace the sensor every 500 cycles.

Remove the short sensor as described in the previous chapter, to reach the situation shown in 4.3.6.

With a size 24 wrench remove the locking nut holding the main nut (Fig.8).

Unscrew the main nut with a 24, hold the sensor cover a size 13 wrench (Fig. 8, fig.9).

With a rubber hammer hit the sensor cover, to put it inside the cell (Fig.10).

Clean the contact surfaces of the sensor cover and cell well.

Use red silicon (temperature resistance till 250°C) to seal the connections between the sensor cover and the cell.

Replace the sensor cover with a new one, mount all parts in place, be sure that all the nuts are tightened well.

Fig. 7

Fig. 8

Fig.9

Fig. 10

WARNING! The sensor cover has influence on the temperature measurement. For this reason use only original spare parts.

4.3.8. Replacing the long temperature sensor.

The long temperature sensor is a part of the protection circuit. For this reason use only original spare parts, to avoid malfunction or a risk situation during use of machine.

- Release the plastic cover, and remove the sensor.
- Disconnect the wires connected to the temperature display.
- Install a new sensor. The polarity is important, so take care about the connection. Watch carefully how the previous sensor was connected, and connect the new one in the same way.
- Tighten the plastic cover.

4.3.9. Replacing the long temperature sensor cover.

Every 1000 cycles replace the cover of the long temperature sensor, following the procedure described below.

- Remove the long temperature sensor, as described in the previous chapter.
- Remove the Teflon cover.
- Unscrew nut number 10 with a size 17 wrench while holding the sensor cover with a size 8 wrench (Fig.).
- With a rubber hammer hit the sensor cover, to put it inside the cell (Fig.).
- Clean the teflon isolation and contact surface of the cover.
- Place a new cover, and put all parts in place. Tighten all the nuts well.

Fig. 11

Fig.12

The long temperature sensor is a part of the protection circuit. For this reason use only original spare parts, to avoid malfunction or a risk situation during use of machine.

4.3.10. Changing the lid gasket.

Every 1000 cycles it is advised to change the gasket of the lid (MC00276), following this procedure.

- Remove the old gasket from mounting ditch.
- Insert the new gasket, begin from point A (Fig.).
- Placing the new gasket in position You can use a screwdriver to push the mounting fin inside the ditch. Be careful not to puncture the gasket Attention the gasket must be cut.
- When the gasket is placed in all the length of the ditch, cut off a fin from the part which is free, and insert the free part inside the fixed side of the gasket.



Fig.13


4.3.11. Replacing the carbon and absolute filter.

Carbon filters, which are placed inside the filter group should be changed every 6 months. The absolute filter should be changed every year.

To replace the carbon filter follow the procedure described below:

To replace the absolute filter:

WARNING! Filters are part of the air treatment circuit of the machine. Use only original or authorized parts.

4.4. At the end of the cycle

At the end of every cycle you must clean the machine.

Using a thoroughly clean wire brush clean:

- the mesh under the lid (close the port that leads to the condensation column);
- the edges of the discharge port and hatch.

The machine will now be ready to start another cycle.

4.5. End of the work session

Cleaning the machine:

Using a wire brush thoroughly clean:

- The grid under the lid (close the port that leads to the condensation column);

- the edges of the discharge port and hatch.

The machine will now be ready to start another cycle.

Leaving the machine unattended:

At the end of the working day it is advisable to leave the machine with the lid and handle locked in place. After having closed them, it is therefore necessary to move the main power switch to the "0" position (the clicks of the locking devices will be heard). To turn the machine off, move the main switch to the "0" position and use the key to lock the door of the control panel. The key must be kept by the operator.

Everyday at the end of shift the operator must carry out simple operations that serve to keep the machine clean and efficient. Following are the operations that must be carried out at the end of each shift:

- Clean the metal net situated on the lid with a vacuum cleaner, if necessary rub the net with a brush and remove the dust on it;
- Clean the bottom of the cell controlling the wear of the blades;
- Check the inside of the cell making sure the screws are tight;
- Clean dust from the nozzles that are located within the cell;
- Disconnect the tube that leads vapours from machine group filters and empty it from condensation that could have formed during operation.

Circular Filter

4.6. Maintenance for long inactivity

If you are planning on not using the machine for long periods of time proceed as follows:

Maintenance

- Switch off general power on the electrical board.
- Thoroughly clean the cell and oil with rust proof oil the bottom of the cell and the rotor blades.
- Clean the filter located on lid with a metal brush and wash the filter with water and leave it in the open air.
- Activate the pump until the circuit is completely empty.
- Empty water present in the tube that connects machine to filter cabinet.
- Protect the machine and all its components from dust.


Finally, in order to keep the machine in perfect working efficiency, even during long periods of non use, we advise that you contact specialised personnel.

5. PROBLEM SOLVING

5.1. Rotor blocked because of overload

If during the cycle the rotor stops because of an overload, proceed as described

in paragraph 3.5 . Now push the  key (The LED will light up and the motor will turn in clockwise) and restart the cycle.

Remember that during the normal functioning it is preferable to maintain the motor in counterclockwise direction (LED on the  key OFF). Use position clockwise direction only to superate conditions blocking the rotar .

There is the possibility with certain malfunctions that the cycle cannot re-start . If it is not possible to re-start and complete the sterilization cycle , the material and the equipment can present biological risks .

Place the machine in "Manual" mode. Open the cell's lid and the discharge hatch, empty the machine into a few bags by pushing the waste through hatch. Load the bags in a new cycle to be reworked.

In the event that a malfunction occurs right after the start of the cycle and the waste has not been sufficiently shredded to be able to extract from the bottom hatch, proceed to extract it by hand from the top after protecting oneself with a smock, heavy-duty gloves and a mask against particulate matter (follow the hygiene regulations closely). At the end of the operation the smock must be washed immediately in water and disinfectant and the other safety equipment must be bagged and treated like hospital waste.

Remember that , by loading big bags containing material that is particularly resistant like, for example sheets, scrubs, rubber and plastic tubes, and if at the same time the blades are not very sharp, it is possible that the resistant torque vastly exceeds that available from the motor. In that event , in fact, the rotor

instead of penetrating the material drags the bags forward without shredding them resulting in overload.

5.2. The printer does not print and the cycle does non finish

If during the cycle : the printer does not print any values , the temperature shown on the display of the operator's panel doesn't rise and at the same time the temperature display on the machine's console shows that the temperature is rising ; the temperature probe is broken. Replace both the temperature sensor and the sensor cover as described in the maintenance chapter.

If the printer still does not print the temperature values after the replacement the cause may be the interface between the sensor and the PLC . In this case call specialized assistance.

5.3. The treated material is too hot when discharged

The temperature allowed for the discharge of the material is 95°C.

If the material is discharged at a temperature of over 120°C control the system introducing water.

Especially the water faucet , the injection nozzles inside the cell and the solenoid valves EV0 and EV1 .

Check the state of wear of the short temperature sensor which commands the PLC.

If after these controls the material is discharged at a temperature above 100°C contact technical assistance.

Attention ! Keep in mind that material that is too hot involves high risks of spontaneous combustion even several hours after being discharged and bagged.

5.4. Strange noise of metal banging inside the cell

The machine is not designed to treat large pieces of metal.

If there are hard, compact metal parts with a weight of over 50 grams, the machine will produce a strange noise caused by the banging of metal pieces against the rotor and the sides of the cell.

In these cases, to reduce as much as possible, the possibility of damaging the rotating and fixed blades, it is necessary to slow down the speed of the rotor ending the cycle in "Manual" mode in low speed.

To do so it is necessary to shut down and restart the cycle in "Manual" mode at low speeds. Complete the treatment by performing all operations (heating, cooling, exhaust) as described in section 3.6, but do not go into high speed. Of course, the duration of the cycle will be longer. At the end of processing open the cell and ensure that there are no more metal masses on the bottom of the cell.

5.5. The machine heats up less than usual

Note that if the cell has been filled and the rotor turns, it is impossible that the machine does not produce heat.

If the temperature remains practically constant and it seems that the warming is weak, check the status of the temperature sensor and of the blades; if the blades are very worn their warming effect is greatly reduced.

It is also possible that the blades are a little worn but the material is packed in bags so light that they "float" and do not descend into contact with the rotor. In this case the material in contact with the rotor is very hot but the average temperature inside the cell remains fairly low. Normally, in these cases, before the end of the cycle the material above reaches the rotor and the temperature rises very quickly. If this phenomenon occurs frequently, the lighter packages must be loaded into the cell first, followed by the heavier ones

However, if these were not the cause, then the fact is due to:

- load of waste is too light

Problem so

- moisture content of the waste is too high
- Remember that the machine is a solid waste sterilizer , so the weight of the material must not be less than 10kg and the percentage of humidity should not exceed 40-45%

5.6. Steam escapes from the cell

The output conduit of the vapors from the cell must remain free so that the suction due to the fan of the filters group maintains the cell in slight depression. In this way, the vapors can not get out of the cell by any other way except through the filter group.

If during the cycle the depression fails, a leak of vapours from the seal of the cell may occur.

It is possible that excessive dust has clogged the circular filter on the lid; in this case it is necessary to open the cell and thoroughly clean the filter and then check that the vapor suction hole is not blocked.

**Vapors suction
hole**

If even after this intervention the steam continues to exit from the cell, put the machine in "Manual" mode and turn on the vacuum cleaner filters in the group, cut a strip of paper 4 cm wide and 20 long and bring it close to the vapour suction hole holding it tightly ; if the paper is not sucked into the tube it is sure that the columns of lowering temperatures are too clogged with limestone. In

this case the intervention of specialized technicians is required , please contact support.

5.7. Table of the main problems, their causes and solutions.

PROBLEM	PROBABLE CAUSES	SOLUTION
Loud squeal coming from the transmission group	Belt not stretched	Tighten the belts and check their condition, change the whole group, if necessary (MAINTENANCE SPECIALIST)
Abnormal vibrations of the machine running in jolts	The rotating blades are not sharp enough or are too rounded on the outside toward the wall	Sharpen or replace the blades, especially if very rounded
	Excessive load of waste	Complete the cycle. In the future, stick to a lighter load
The time required to cool the treated material is excessively long	The nozzles that send water to the cell are blocked	Clean the nozzles
	EV0 and EV1 valves damaged	Check valves functionality (MAINTENANCE SPECIALIST)
	Inconsistency or lack of water coming from the water supply	Check that there is pressure in the water supply and that there are no shutoff valves closed
The treated material has a bad smell such as organic type	Both carbon filters and absolute filter need to be changed	Replace both filters inside the filters group (MAINTENANCE SPECIALIST).
The treated material has a strong smell of toast	The material was kept at too high temperature (approx. 200 ° C) due to dysfunction in the system to measure the temperature	Check the status of the temperature sensor
Loss of vapors from the cell	Clogging of the duct of the vapors that connects the cell to the absorption columns	Clean the duct.
	Filter group fan not working	Restore the fan functionality (MAINTENANCE SPECIALIST)

PROBLEM	PROBABLE CAUSES	SOLUTION
	Both carbon filters and absolute filter are clogged	Replace both filters inside the filters group (MAINTENANCE SPECIALIST).
	The lid gasket is damaged.	Replace the gasket
Loss of vapors from the cell and loss of water under the machine	The water level in the first absorption column has exceeded the overflow level and reached the height of the duct of the vapors; in these situations, there is no suction	Check that there are no blockages in the column. Remove sediment on the bottom of the tank; decrease, if excessive, the flow of water sent to the column (MAINTENANCE SPECIALIST)
	The circular filter is clogged	Clean the circular filter
	The conduit of the vapors that connects the cell to the absorption column is clogged	Open the lid and remove the deposits in the vapors conduct
	The filter fan is not working properly or not working at all, clogged filters	Check the filter fan and the connection with the absorption column; check the air filter and replace it if necessary (MAINTENANCE SPECIALIST).
It seems that the machine heats less than usual	To much water in the waste load or waste too light	Prepare the load more evenly
	Loss of water in the cell for failure of the EV1 and EVO valves.	Remove and clean the valves (MAINTENANCE SPECIALIST).
	Both carbon filters and absolute filter are clogged	Replace both filters inside the filters group
	The electric heaters that heat the walls of the loading mouth do not work.	Check the fuses and replace the faulty resistors.