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
The image shows a 'Login' dialog box with a close button (X) in the top right corner. It contains two input fields: 'User:' and 'Password:'. At the bottom, there are two buttons: 'Cancel' and 'OK'.

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 A screenshot of a 'Login' dialog box. It has a title bar with the word 'Login' and a close button (X). The main area contains two input fields: 'User:' and 'Password:'. At the bottom, there are two buttons: 'Cancel' and 'OK'.

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


Fig. 55

| CELL GROUP, GRNW050-04-02 | |
|----------------------------------|---------------------------------------|
| CODE | DESCRIPTION |
| A500279 | Piston 0.20 C.15 Single Effect |
| FC00199 | Washer Zn 10.5x30x5 |
| FX00229 | Washer 24,5x39x3 Stainless |
| FX00261 | Screw TCEI Inox M8x25 |
| FX00450 | Nut M15 Short Stainless |
| FX00452 | Nut M24 Short Stainless |
| FX00455 | Washer 17x23x3 Stainless |
| FX00457 | M10 NUT A2 SHORT |
| M000289 | Gas Spring Mod. 198.390.35 |
| M000391 | Seal Edge Cell |
| ME00456 | P1100 Probe |
| M100517 | Vassel Sleeve Sp.35 |
| M000739 | Cover Short Sensor 35 |
| MC00751 | Cover Long Sensor Vassel 35 |
| MC00752 | Protection measuring hole Vassel 35mm |
| MC01588A | Protection Vassel NW50 |
| MC01588B | Protection Vassel NW50, Sensor Holden |
| MC01612 | V-Ring Seal V-110A |
| MC01658 | Vessel ProteEtion 2nd lever, Hardox |
| MC01659 | Discharge Mouth Protection NW50 |
| MC01663A | Bottom Vassel Protection Hardox NW50 |
| MF00742 | Sensor Short Temperature |
| SANW050-02-02 | Lid NW50 |

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S.A. ACQUA-WATER SUBASSEMBLY, SANW050-24-00

Fig. 60

| SANW050-24-00 | |
|---------------|------------------------|
| CODE | DESCRIPTION |
| MF00843 | Heat Exchanger NW50 |
| MF00842 | Exchanger below part |
| FX00287 | Washer Inox 6,5x18x1,5 |
| MF00844 | Exchanger top part |
| FX00285 | Nut M6 Stainless |
| FX00220 | Grower Washer Zn M6 |

| | |
|---------|-----------------------|
| MC00547 | Spacer Bush |
| FX00295 | Screw TCEI Inox M6x20 |
| MC01741 | OR CORDA 250,2 x 5,34 |

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| | | |
|-----------|--|-----------|
| 3.2.7. | Language selection | 41 |
| 3.2.8. | Set up of cycle parameters in the sub menu | 41 |
| 3.2.9. | Inverter Amper setup..... | 43 |
| 3.2.10. | Hertz Settings | 43 |
| 3.2.11. | Inverter Scale | 44 |
| 3.2.12. | Temperatures of the cycle | 44 |
| 3.2.13. | Hospital data and serial number | 46 |
| 3.2.14. | Bypass inspection door | 46 |
| 3.2.15. | Timer of the cycle | 47 |
| 3.2.16. | Setup sensor data | 50 |
| 3.2.17. | Manual Function | 51 |
| 3.3. | Switching on and switching off..... | 51 |
| 3.4. | Automatic cycle..... | 52 |
| 3.5. | Manual machine controls..... | 52 |
| 3.6. | Service messages..... | 54 |
| 3.7. | Alarms..... | 54 |
| 3.8. | Filling of sterilization vessel | 59 |
| 3.9. | Carry out a pseudo-cycle in 'Manual' | 59 |
| 3.10. | Downloading material from the cell in MANUAL mode..... | 61 |
| 3.11. | Verifying the Effectiveness | 61 |
| 3.12. | Integrated Waste Collection system..... | 62 |
| 4. | Maintenance | 64 |
| 4.1. | Informative notes..... | 65 |
| 4.2. | General norms..... | 65 |
| 4.3. | Replacing consumable material | 66 |
| 4.3.1. | List and duration of consumable parts | 66 |
| 4.3.2. | Changing printer paper | 77 |
| 4.3.3. | Vertical blade changing. | 77 |
| 4.3.4. | Blade support dismanteling..... | 78 |
| 4.3.5. | Changing blades..... | 79 |
| 4.3.6. | Replacing the short temperature sensor. | 79 |
| 4.3.7. | Replacing the short sensor cover. | 80 |
| 4.3.8. | Replacing the long temperature sensor..... | 81 |
| 4.3.9. | Replacing the long temperature sensor cover..... | 81 |
| 4.3.10. | Changing the lid gasket. | 83 |
| 4.3.11. | Replacing the carbon and absolute filters. | 84 |
| 4.4. | At the end of the cycle | 88 |
| 4.5. | End of the work session..... | 89 |
| 4.6. | Maintenance for long inactivity..... | 89 |
| 5. | Trouble shooting | 91 |
| 5.1. | Rotor blocked because of overload | 92 |
| 5.2. | The printer does not print and the cycle does non finish | 92 |

5.3. The treated material is too hot when discharged 92

5.4. Strange noise of metal banging inside the cell 93

5.5. The machine heats up less than usual 93

5.6. Steam escapes from the cell 94


5.7. Table of the main problems, their causes and solutions..... 94


1. General information

1.1. Identification data of the manufacturer

Newster System s.r.l.
Via Pascoli, 26/28, 47852 Cerasolo, Coriano (RN) - Italy
Tel. +39.0541.759160
Fax. +39.0541.759163
web site : <http://www.newstergroup.com>


1.2. Identification of main parts

Vapour to the heat exchanger

Fig. 2 - Cell Group

Heat

Vapour inlet

Water inlet

Vapour outlet

Water outlet

Fig. 3 - Demister Group

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
Fig. 5 - Filter Group

Fig. 6 - Discharge Vessel

Fig. 7 - Inverter Group

Fig. 8 - Electrical Board

Fig. 9 - Oil Pump

Fig. 10 Pump and water tank (optional)

Fig. 11 Radiator (optional)

1.3. Identification plate

For any communication to the manufacturer or the distributor, it is necessary to refer all the information on the identification plate.

Fig. 12

The identification plate, including the machine serial number, is located on the top of the machine, in a visible position.

1.4. Use of the machine

The purpose of the NEWSTER®NW30 is to treat potentially infectious solid healthcare waste in order to obtain:

- Sterilization
- Physical modification

- Dehydration
- Volume reduction
- Weight reduction

The machine and its components must be used only for the treatment of solid healthcare waste, containing liquids, potentially contaminated by pathogenous microorganisms, deriving from hospitals, laboratories, medical research centers, dialysis centres and infectious disease hospital wards. Kind of health wastes treated by Newster and their classification:

| European code | Kind of waste | Classification |
|--|--|-------------------------------|
| C.E.R. 180103 or 180202 Wastes coming from research activity and bacteriological diagnostic | Sanitary towels, napkin and nappy Cotton flock for papetest and colposcopy Ocular sticks not sterile TNT ophthalmic sticks Cannulaes and drain pipes Catheters (vesical, venous, arterial, pleural drains, pipe-unions, probes Extrabody circulation circuits Expired blood bags, blood bags, blood Defluxors Contaminate phleboclysis Dialysis filters, exhausts air filters (without chemical risk) Surgery gloves Disposables : vials, pipettes, tubes, protective dresses , mask, glasses, lenght of material, sheets, foot wears, white overalls Materials for medication (bandages, tampons, wax plasters) Bags (transfusionals, nutritionals, urine) Infusion set Rectal and gastric probes Probes (nose, broncoaspiration, oxygen therapy) Auricular speculum Vaginal speculum Plaster cast or bandage Teeth and small part of body not recognizable Bedstead for experimental animals Empty containers Wastes coming from dentists Wastes from restaurant (infective only) Petri plates, and all contaminated disposables | Dangerous with infective risk |
| Cutting wastes C.E.R. 180103 or 180202 | Needles, syringes, knives, venflon, razors and bistoury disposables | Dangerous with infective risk |
| Cutting wastes not used C.E.R. 180101 or 180201 | Needles, syringes, knives and bistoury | Special |

| European code | Kind of waste | Classification |
|---|---|--|
| Anatomic wastes C.E.R. 180103 or 18202 | Tissues, organs and anatomic parts not recognizable Animals for experiments | Dangerous with infective risk |
| Empty containers | Empty containers from drugs, veterinary drugs, disinfectants, food, drinks and infusion solutions | Special/assimilable to civil wastes |

The machine is not suitable for treating:

- Gas bottles or similar
- Chemical products
- Inflammable material
- Explosive material
- Radioactive material
- Stones
- Bulky or hard pieces larger than 100 gr
- Any material similar to the above listed.

The efficiency of the sterilization system is proved according to the technical norm UNI 10384.

1.5. Installation of the machine

The machine and its accessories could be installed only by the manufacturer or authorized people.

Attention! The electrical connection must comply with the rules of the country where the machine is installed. For this reason, the cable IC00528 must be considered part of the link but can not be the only part of connection to the electricity network. There must be also a protection, a connection to the earth and insulation in accordance with applicable regulations.

1.6. Technical specification

| | |
|-----------------------------------|---|
| Nominal waste disposal potential: | 70-80 kg/hour of waste depending on the humidity and composition of the waste 700 liter/hour |
| Final volume of treated waste: | 20-25% of initial volume |
| Final weight of treated waste: | 70-75% of initial weight |
| Volume of sterilization cell | Approximately 400 liters |
| Overall control system: | Programmable logic controller (PLC) |
| Water consumption: | Approximately 50 liters per day with the water recycling system |
| Electrical power requirements: | Up to 80 kW. |

| | |
|---|---|
| Connections: Electrical Water supply Water discharge | 50/60 hz, 400 Volt, 3F+N+ ground Connection to water grid with a ¾" tap Connection to sewage system by a 50 mm pipe |
| Dimensions: Sterilizer Filter storage Electrical board Inverter group | 250 x 140 x 180 cm 80 x 50 x 110 cm 80 x 45 x 110 cm 80 x 50 x 180 cm |
| Weights: Sterilizer Filter storage Electrical board Inverter group Accessories Total weight: | 2200 kg 70 kg 100 kg 150 kg 100 kg 2200 kg |

Emissions:

The air emissions discharged from the machine conform with the environmental protection requirements, as witnessed by the certification located at the headquarters of the manufacturer.

Before being discharged, the air is treated with active carbon filter and then with an absolute filter. The absolute filter has an efficiency DOP >99,99%.

The water that is discharged in the sewer conforms to the requirements, as witnessed by the certification located at the headquarters of the manufacturer.

The noise level is cut down by polyurethane shields.

The machine does not produce any waste other than the waste treated.

1.7. Disuse of the machine

When the machine is not in use, all the machines components must be disposed of conforming to the rules and regulations of the country where the machine is installed. The recyclable parts must be separated and disposed of at the proper disposing facility. Electrical components must be separated from the rest of the material.

1.8. Warranty

For warranty conditions see the attached form that must be filled and sent to Newster System Srl.

2. Safety notes

2.1. Preliminary observation

Access to the rooms assigned for the use of the machine must be permitted only to personnel authorized by the healthcare management authority.

You are reminded that in the rooms assigned for use of the machine it is forbidden to drink, smoke or eat.

For all operations required, the use of protective gloves is compulsory.

The various risk situations that might affect the machine have been analyzed, verifying the precautions and technical solutions required to eliminate or reduce them in the most appropriate way by applying adequate measures of protection.

2.2. Proper use of the sterilizer "Newster® NW30"

The device and its components may only be used for the treatment of solid potentially infectious hospital waste as listed in paragraph 1.3.

Do not use the equipment to treat gas containers, chemical, flammable, explosive or radioactive materials. Hard and compact masses (eg. compact metal bodies with mass greater than 100 grams, stones, etc..) can cause damages to the rotating and fixed blades.

The apparatus is not suitable for the treatment of liquids or waste having less than 60% of dry matter. In these cases, the waste may be treated after being mixed with material from a previous treatment.

The apparatus can not be loaded only with materials with a melting point of less than 150°C (as PVC). In this case the waste can be treated by mixing with another material which has a melting point higher than 150°C. In any case, the average weight of materials with a melting point less than this value should not exceed 60% of the total weight.

The equipment is designed to operate at atmospheric pressure (in order to enable the cell to be in depression). Do not cover, for any reason, the exhaust outcoming vapours, otherwise the system may go slightly under pressure.

2.3. Qualified Personnel

Operating on this equipment is allowed only to trained personnel, namely:

- operational staff who has received proper training to work with the equipment and who knows the contents of this manual and of the annexed documentation,
- service and maintenance personnel shall be trained to repair automation equipments and authorized to connect and disconnect from electrical power, open, ground and work on electrical circuits, components and systems according to what is established by safety regulations.

Do not leave the equipment available to persons who are not properly trained. The non-compliance of the information contained in the manual or affixed on the equipment itself can cause serious damage to people and property. After opening the panels of the machines and / or of the electrical board panels, certain parts of the machine become available, they can

have dangerous voltages, which could result in death or serious injury to persons and property if appropriate precautions wouldn't be taken.

2.4. Working area


Fig. 13

2.5. Protection against accidents

The device is designed with the following protection against accidents. They should not be removed for any reason:

- electrical equipment follow the standards EN 60204-1
- the auxiliary circuits are made in low voltage (24 V),
- sterilization of the cell operates under a pressure lower than atmospheric .
- the main electrical disconnecting-switch and electrical board door are interlocked,
- an electrical interlock prevents the starting of the rotor if the lid is not closed properly,
- an electrical interlock prevents the opening of the lid if the rotor is in motion,
- to close the lid, the operator must first unlock a safety piston,
- the system provides a self-check to verify the closure of the discharging door,
- in case of lack of tension, the lid of the cell and the discharging door remain tightly closed,

- in case of failure the machine stops in the best conditions of hygiene and safety allowed by the situation.

2.6. Personal Protective equipment

The operator during normal operation of the machine and during the inspection and routine cleaning must use the following personal protective equipment:

- Mask in accordance to EN 149

Fig. 14

- Glass in accordance to the EN166 e EN170

Fig. 15

- Disposable gloves

Fig. 16

- In the case of ordinary or extraordinary maintenance is recommended to use anticut protective glove in accordance to EN 388

Fig. 17

2.7. Labels and safety symbols

Labels bearing symbols and/or safety messages are applied on the sterilizing machine. In this section we show their exact location and the kind of hazard.

Attention. Read the user manual carefully before using the machine or doing maintenance.

Fig. 18

Danger. Moving parts. Danger of crushing arms during operation of the main motor. Never put your hands under the sterilizing machine while the main engine is working.

Fig. 19

 **Danger.** Risk of crushing upper limbs during the inspection and cleaning of the discharge mouth. Never put your hands inside the inspection hole of the discharge during the closing of the mouth. Close the inspection hole by acting on the outlet with both hands. To clean the drain use the tool provided with the machine (see figure 11).

Fig. 20

Danger. Hot surface. Danger of burns of the upper limbs during the operation of the machine. Do not hold the lid by pulling on the metal parts, open and close the lid acting only on the plastic handles.

Fig. 21

Attention. Make sure that all the labels are always readable. Wipe them with a cloth with soap and water. Do not use solvents. Replace the damaged plates receiving them from the manufacturer. If a plate is located on a piece that must be replaced ensure that the new one is of the same type and is located in the same position.

Danger. Risk of electrical hazard due to presence of dangerous residual tension with main switch open for 5 minutes. Absolutely do not work on the electrical equipment before the specified time.

Fig. 22

2.8. Residual risks

The careful analysis carried out by the manufacturer and stored in the file, made it possible to eliminate most of the risks associated with the conditions of use of the steriliser. The manufacturer recommends to strictly follow the instructions in this manual and safety regulations, including the use of protection devices required.

The residual risks can be:

Electrical hazards during maintenance under tension

Following the DLgs 17/2010 only trained personnel is allowed to operate under tension and in accordance with the double insulation to ground.

Risk of improper use in the presence of explosive and flammable substances

Follow the instructions in paragraph 2.2.

Risk of operating in a wet environment

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Risk of crushing of the upper limbs during the cleaning of the discharge

The machine has a label indicating the security risk. It is possible to avoid this risk by using the special tool as shown in section 2.6.

Risk due to moving parts

The machine has a label indicating the security risk. It is possible to avoid this risk by acting as specified in paragraph 2.6.

Risk of burns of the upper limbs

The machine has a label indicating the security risk. It is possible to avoid this risk by using the special plastic handles as indicated in paragraph 2.6.

Risk of misinterpretation of the pictograms

The symbols are installed in accordance with the rules for graphic symbols to be used.

The user must immediately replace the safety signs that may become unreadable due to wear or damage.

Attention. It is absolutely forbidden to remove the safety labels on the machine. The manufacturer declines all responsibility in the event of breach of this prohibition.

The machine described in this manual applies to the following guidelines: 2006/42/EC, 2004/42/EC and 2006/95/EC.

2.9. Authorized personnel and responsibility

To prevent injury, this apparatus must be managed, by the person responsible, one operator at a time.

When the equipment is turned on, the support staff and maintenance must remember not to put their hands under the edges of the lid of the cell and must perform operations only after obtaining the consent of the responsible conductor at that time.

2.10. Hygiene norms

Always remember that the handling of hospital waste involves the potential of incidental contact with pathogenic microorganisms. So: at work wear the protection clothes (to be replaced at least weekly and washed with disinfectant) and use disposable gloves when handling bags containing waste. If the wastes are originated from infectious wards, increase caution during these operations, by also wearing a mask, a pair of glasses and respiratory protection.

Keep the working space well ventilated (at least 10 air changes per hour).
Do not smoke, keep or consume food in the workplace.

Operators and the technical assistance should not touch their face, especially the mouth, nose and eyes before washing hands with disinfectant soap.

Always wash your hands with disinfectant soap after work and before eating, drinking, smoking (in a separate room).

Even though before the end of each cycle the equipment sterilizes itself, for caution before performing maintenance on parts which can come into contact with the waste or vapors, sprinkle abundantly with disinfectant (disinfectants used: based on chlorine, iodine, hydrogen peroxide and the like), close the lid and turn on the heating elements, wait 15 minutes for the disinfectant to act, turn on the vacuum cleaner filters and after 5 minutes turn everything off and then begin to maintenance work. At the end of maintenance work the protective coat needs to be changed.

In case you have an accidental spill of liquid on the floor, proceed to the collection and cleaning with disinfectant.

At the end of the work, wash the floor with water and disinfectant.

2.11. Protection from electric shock

The equipment must be connected to a good grip on the ground before being turned on.

The operation and maintenance of electrical components and automation equipment may be carried out only by trained personnel.

The power cables to the board are always in tension. For eventual maintenance to the input clamps of the electric board the operator must remember to interrupt the power supply to the electric board

Never throw water on live electrical components.

Always turn off the machine at the end of the working day and in any case before maintenance you must make sure that the electrical circuit is completely de-energised by turning the disconnecter, located on the left of the front panel, horizontally to 'O' (Or 'OFF').

2.12. Emergency stop

In the case of urgent need to stop the device, press the button '**EMERGENCY STOP**' obtaining the immediate arrest of all the components, for the fall of the auxiliary circuits.

The buttons "EMERGENCY STOP" are two and they are: in the lower right corner on the console of the machine (Figure 19) and on the control panel as shown in (Figure 20).

Fig. 23

Fig. 24

2.13. Electric welding on the equipment

Electric welding on the structure or any part of the equipment can cause serious damage to electronic components. Request for service and disconnect the PC from the electrical panel.

2.14. Operations on components sensitive to static electricity

Do not allow unqualified personnel to do maintenance to electronic components. The basic rule to follow is that electronic components should be touched with the hands only if you need to work to be performed on them. Do not touch the components contacts and their conductors.

Touch components only if your body is electrically discharged continuously to the ground by means of a wrist strip conductor or special shoes or shoes with antistatic conductive strips to the ground.

Before touching an electronic module, be sure to not be electrostatically charged. The easiest way is to touch something grounded (such as a water pipe) before touching the module.

2.15. Preventing fires

Keep a dry powder fire extinguisher on hand. Water extinguishers are not suitable..

Never spray water on electrical components.

To avoid risk of combustions, before leaving the workplace unattended, ensure that the treated material is not hot over 100 ° C.

If the material is discharged at temperatures above 130 ° C, it must cool with water and it must absolutely not be close in the collection bag until its temperature is not less than 80 ° C. A material treated with too high temperature could cause an a fire.

If the material being processed has a strong smell of roasted, it means that it has been subjected to excessive temperatures (over 180 ° C) due to malfunction. The most likely causes are the lack of water during the cooling phase, or anomalies to the system of temperature measurement. Control the transport system of the water and the chain of temperature regulation.

2.16. "Logical" equipment

The normal operation of the equipment is in "Automatic" mode via a logic and a sequence of phases and interlocks that follow one after another. This logic is factory set and is contained in the EPROM of the CPU of the PLC. The times can be changed only by specialized programmers (if authorized and provided the password) and are held in RAM memory by the lithium battery. If you remove the battery to restart it will load the default time data stored in the EPROM from the manufacturer.

Changes to the electronic logic involves a risk of malfunction and may only be performed by experienced programmers in automation.

2.17. Locked rotor stress

No particular caution must be taken during the normal operation of the machine because its interior is accessible only when the system has passed through a state of sterility. However, in case of locking without possibility of restarting before the completion of the sterilization cycle, it is necessary to use some caution.

Wear a lab coat, heavy-duty gloves (eg. Leather) and a mask with filter (follow the hygiene rules in this manual) when you open the cover always be protected by gloves, glasses and mask and empty the machine by pushing the material through the exhaust port in a few bags that can later be reworked into a new cycle.

If the failure occurred immediately after the start of the cycle, and if the material is not yet sufficiently shredded to be extracted through the discharge port, it is unfortunately necessary to proceed manually from above. Wear a lab coat, heavy-duty gloves (eg. Leather) and a mask with filter (follow the hygiene rules in this manual).

2.18. Operating precautions

All control devices are clearly visible, and the operating cycle has been planned to simplify and limit the activity of the operator to the turning on of the machine, the loading of the waste into the sterilization vessel and the starting of the treatment cycle.

To avoid all risks, operation in automatic mode (cycle selector set to automatic) for the entire cycle provides a suitable guarantee. However, it is also possible to select manual mode for the cycle to check operation without load or to solve problems that arise during automatic operation.

The risk of access to the operating zone of the machine (vessel with rotor and blades) has been eliminated by:

- system that locks out start of cycle if lid is not closed (with circuit microswitch);
- system that locks out start of cycle if discharge port hatch is not closed (with circuit microswitch);
- system that locks lid in place during the cycle of the machine;
- system that blocks the opening of the discharge hatch during the operation of the machine.

The potential risk of breakage of the rotating blades or parts of the rotor support system (rotating flange and blade holder) are analyzed under normal operating conditions.

The components of the rotor are sized according to the stress envisaged for the treatment by pulverization of materials currently used in hospitals (in particular needles, scalpels, sanitary boxes, metallic probes, glass, as well as plastics, cellulose products, organic wastes, etc). The rotor is therefore capable of resisting the impacts and stress foreseen for its operation.

However, to avoid damage or injuries from the rotating blades, due to the centrifugal force to which they are subjected, in case of slackening or the escape of the blade from the blade holder on the rotor, periodic checks on the tightness of the blades to the blade holder of the rotor are advised.

Using a suitable torque wrench check the tightness of the six bolts holding the blades to the rotor, applying a tightening torque of 9 kgm to them, we advise you check tightness at least every 1000 operating cycles.

Every time the blades are removed the fixing bolts must be replaced. The replacement bolts must have a breaking strain greater than 80 kg/mm².

The machine is equipped with a button for normal stopping and a suitable emergency stop button with a manual reset system to return the machine to operational status.

The following possible malfunctions of the machine during operation were considered:

- In case of a power failure, the lid and the discharge valves for materials remain locked in place to prevent their opening or the emission of materials from the vessel when the machine cycle starts again
- In case of partial malfunctions after the start of the cycle (lack of pressure in the water circulation system, partial electrical malfunction, failure to register temperature, etc) the machine is stopped and the problem is warned of on a specific display, so that it can be rectified before the automatic cycle is restarted.
- The control and monitoring system is easy to use, with easily identifiable buttons, simple instructions and alarm messages shown on the display or with indicator lights, with ample explanations in manual.

The electrical system of the machine is compliant with the CEI EN 60204-1:2006 standard, harmonised in the EU zone. The machine must be used only when connected to a ground

cable and system, as envisaged in the installation phase (inside the door of the electrical panel). The power supply line must be protected with a delayed-action circuit breaker set to trip above 100 A. Power supply cables for the machine must have a cross-section of 16 mm² if cable length is shorter than 20 meters, 20 mm² if cable length is between 20 and 50 meters, and 30 mm² if cable length is over 50 meters. Electrical installations must be carried out by a qualified expert in electrical systems.

At the end of daily operating sessions the electrical system of the machine must be turned off with the main switch located on the electrical panel (position "0"), and this operation must normally be carried out by the user of the machine.

To avoid the risk of spontaneous ignition, at the end of daily operating sessions it is essential to check that any treated materials remaining in the vessel are at a temperature of less than 100°C.

When handling healthcare wastes, which create the hazard of contact with pathogenic micro-organisms, the operator must wear protective overalls and gloves to handle the bags containing wastes. The overalls must be made of a fabric resistant to acids and alkalis, and must be changed at least once weekly and sent for hot washing in chlorine bleach. Gloves (compliant with EN 374) must be of the disposable type. If wastes from isolation wards are to be treated, when the vessel lid must be opened or there are wastes or residues inside the vessel during internal vessel cleaning operations, all handling and movement operations must be carried out wearing a face mask (compliant with EN 140) to protect the respiratory tract, in addition to overalls and gloves.

The room in which the machine is installed must be ventilated with an extractor that performs at least 10 changes per hour of the total air volume. The room must also be equipped with a washbasin, fire extinguisher and running water. The operator must be warned to refrain from touching the face (mouth, nose and eyes) unless the hands have been washed first. At the end of every work session it is always compulsory to wash hands with a disinfectant soap (i.e. with povidone-iodine like Betadine or a similar product). Provide warning signs and forbid smoking and the keeping of food in the room.

For the machine described in this manual these directives can be applied: 2006/42/CE, 2004/42/CE e 2006/95/CE.

2.19. Maintenance precaution

For maintenance operations, and in particular for maintenance operations inside the vessel, the main switch must always be turned off.

To deactivate the electrical system the main switch on the electrical panel must be moved to the "0" position.

Maintenance operations on the electrical system must be carried out by specialized personnel, as the electrical panel is interlocked with the opening lever of the panel, although the power supply terminals still have energy if the electrical panel is opened.

All maintenance operations must be carried out using the personal safety equipment envisaged for this type of work. At the end of maintenance operations, protective overalls must be changed.

A dry-powder extinguisher with a capacity of 6 kg must be kept in the vicinity of the electrical panel of the machine or at some other accessible point in the room where the machine is installed, for use in case of possible emergencies.

2.20. List of required precautions and equipment

Personal safety equipment (to be worn in its entirety during maintenance operations on parts in direct contact with wastes):

- protective overalls resistant to acids and alkalis (to be hot-washed in chlorine bleach once weekly);
- disposable gloves (compliant with EN 374);
- face mask (for cleaning vessel compliant with EN140);
- protective glasses and eyewash.


2.21. Other precautions:

- Washbasin with running water and disinfectant soap; availability of a water can;
- dry-powder extinguisher capacity 6 kg;
- warning signs (including prohibition of smoking and storage of food in room);
- information to operator about hygiene precautions to be taken (for mouth, nose, eyes);
- ventilation of room (10 changes of air volume per hour).

3. Instructions for the user

3.1. Description of process

The NEWSTER®NW30 converter-sterilizer is designed to treat Healthcare Risk Waste (HCRW) at atmospheric pressure and high temperature in a wet environment, as illustrated in the following diagram:

- First stage : loading and start of heating phase
- Second stage : evaporation of liquids
- Third stage : superheating and sterilization
- Fourth stage : cooling
- Fifth stage : unloading

In a closed sterilization vessel a powerful rotor fitted with blades disintegrates, agitates and heats the wastes by impact, friction and resistances.


The temperature of the mass of wastes is measured in real time and with great accuracy by special sensors, while they are being constantly and vigorously agitated.

When the temperature reaches the predetermined level of 150°C, the mass of wastes is automatically sprayed with water so that the temperature decreases.

The treated wastes are cooled down to 95°C. At this point the cycle has been completed, and the product, by now sterile, is automatically unloaded.

To disperse the heat produced by the system, part of the water is continuously replaced by fresh water from the mains supply. Excess water and incondensable gases are discharged into the sewers, with values within legally established limits.

The process lasts about 30 minutes, depending on the ratio of organic matter in the waste, and automatically follows these stages:

1. Wastes are loaded into the sterilization vessel, the lid is closed and the treatment process is started by pressing a button. The engine starts in the first speed, filters start and resistances are turned on.
2. At 60°C the general water electrical valve opens and water starts to flow in the cooling columns. The rotor turns slowly at first, starting to pulverize the materials, and at the same time the temperature starts to rise.
3. The rotor turns faster (going from the first speed to the second), the temperature starts to rise rapidly and the materials are finely pulverized.
4. When 96–100°C are reached, the temperature remains stable until the water present in the wastes has completely evaporated.
5. After the water has evaporated the temperature starts to rise rapidly again, reaching 150°C. After the peak the rotor turns slowly (first speed) and resistances turn off.
6. The mass of wastes is sufficiently moistened by a water spray to cool down to 95°C.
7. The sterilization cycle has now been completed. The vessel is opened, and the product is extracted and collected in the stainless steel integrated waste collector.

3.2. Comands


Fig. 25

3.2.1. The operator panel

The normal operation of the equipment will take place in “automatic” mode by means of a logical linking of phases and blocks that follow each other in an automatic sequence. The logic is pre-loaded in the factory and is contained in the programmable logic controller (PLC) on an EPROM card. This card contains pre-loaded default operating parameters. Only an authorised user with the key software is able to enter into a dialogue with the system and to introduce changes to the parameters.

WARNING! Modifications to the electronic logic carry the risk of disfunction and thus making the guarantee void. In any case, such modifications must only be made by programmers who are experts in automation.

The operating panel permits the control and variation of the timers and temperature values which control the normal course of the automatic cycle. It can also send alarm messages and information on the state of the machine.

The operating panel has a touch screen (figure 26)

Fig. 26

3.2.1. Main menu

By pressing the icon  at any time the main menu will appear.

By pressing the icon  at any time the alarm messages menu will appear as described in the proper chapter.

Fig. 27

You can check the status of the cycle at any time by touching the  icon:


Fig. 28

From the main menu you can enter 4 submenus.

Icon  to display the inner cell temperature (par. 3.2.2)

Icon  to display the system submenu (language, date and time, input/output etc) (par. 3.2.3)

Icon  to display the manual function submenu (par.3.2.16)

Icon  to display the cycle parameters set up in the submenu (par.3.2.7)

Icon 

Icon  to display the bin washing machine.

3.2.2. Display the inner cell temperature

Use this page every time you want to know the temperature inside the cell. Temperatures below 40°C will not be shown.


Fig. 29

Press the icon  to go back to settings submenu.

3.2.3. The submenu settings


Fig. 30

Icon  to display the panel contrast set up in the submenu (par.3.2.4)

Icon  to display the INPUT and OUTPUT panel (Figure 26):


Fig. 31

Icon  to set date and time in the submenu (par. 3.2.5).

LANGUAGE

Icon  to change the language in the submenu (par.3.2.6).

Other function buttons are used only by expert technicians. Refer to the user manual of the PLC KPT700-KP8 Siemens to find out about these functions.

If you have questions regarding the use of the products described in the manual and are unable to find an answer here, please contact your nearest Newster representative.

3.2.4. Set up contrast


Fig. 32

Press the icon of the desired display contrast value to set the contrast.

**CALIBRATE
DISPLAY**

Press the icon  to calibrate the touch screen

3.2.5. Set up time and date



Fig. 33

Touching the blank field you want to fill out this window will appear:

Touch again on the blank field an alphanumeric keyboard appears. Use it to write all the data.

Press the icon  when finished, and to start the clock.

Press the icon  to go back to settings submenu.

3.2.6. Language selection


Fig. 34

Select the language by pressing the language chosen icon; PLC changes the language immediately.

3.2.7. Set up of cycle parameters in the sub menu


Fig. 35

Press the icon  to display the inverter amper set up (par. 3.2.8)

Press the icon  to set up temperatures of the cycle (par.3.2.11)

Press the icon  to set up hospital data and serial number (par.3.2.12)

- Press the icon  to set up inspection door switch (par.3.2.13)
- Press the icon  to set up the timer of the cycle (par.3.2.14)
- Press the icon  to set up the sensor data (par.3.2.15)
- Press the icon  to display the total number of cycles and the number of cycles on the printer (figure 31.)


Fig. 36

- Press the icon  to display the inverter hertz set up (par.3.2.9)
- Press the icon  to display the inverter scale set up (par.3.2.10)

3.2.8. Inverter Amper setup


Fig. 37

Press the icon  to go back to cycle parameters settings submenu.

3.2.9. Hertz Settings



Fig. 38

Touching the blank field you want to fill out this window will appear:

Touch again on the blank field an alphanumeric keyboard appears. Use it to write all the data.

Press the icon  to go back to cycle parameters settings submenu.

3.2.10. Inverter Scale



Fig. 39

Touching the blank field you want to fill out this window will appear:

Touch again on the blank field an alphanumeric keyboard appears. Use it to write all the data.

Press the icon  to go back to cycle parameters settings submenu.

3.2.11. Temperatures of the cycle


Fig. 40

Press the icon  or  to move between the set up temperatures of the cycle submenus pages.

Press the icon  to go back to cycle parameters settings submenu.

In these submenus you can set up the temperatures of the cycle, furthermore these temperatures are printed to describe the course of the cycle. Not all values are possible, inside PLC are configured some limits and every temperature must respect these limits.

It is highly recommended to contact the producer before changing these values to avoid any malfunction during the cycle.

Severe personal injury or substantial property damage can result if proper values are not entered in this data base.

Touching the blank field you want to fill out this window will appear:

Touching the blank field again an alphanumeric keyboard appears. Use it to write all the data. In this data base it is possible to set the following values.

- Opening general water: general water valve opens at this temperature and water starts to run inside heat exchanger column.
- Pump restart: pump with deodorizer (where it is present) starts again at this temperature.

- Sterilization temperature: it is the maximum temperature reached during the cycle
- Closing water in cell: at this temperature the discharging process starts during the cycle

3.2.12. Hospital data and serial number



Fig. 41

Press the icon  to go back to cycle parameters settings submenu. In these submenus you can set up all the hospital data and the serial number of the sterilizer. This data is printed on the reports. The hospital data must stay on 2 lines of 20 characters each; meanwhile the serial number can have only 3 numbers. Touching the blank field you want to fill out this window will appear:

Touch the blank field again and an alphanumeric keyboard appears. Use it to write all the data.

3.2.13. Bypass inspection door



Fig. 42

Touching the slide  icon this window will appear:

Touch the blank field again and an alphanumeric keyboard appears. Use it to write username and password. Now it is possible to activate or deactivate the safe alarm for inspection door and for the recirculating system.

We warn that is possible to bypass this alarm only for maintance by authorized operator.

3.2.14. Timer of the cycle

Fig. 43

Fig. 44

Press the icon  or  to move between the set up temperatures of the cycle submenus pages.


Press the icon  to go back to cycle parameters settings submenu.

In these submenus you can set up all the timers of the cycle. Not all values are possible, inside PLC some limits are configured and every timer must respect these limits.

It is highly recommended to contact the producer before changing these values to avoid any malfunction during the cycle. Severe personal injury or substantial property damage can result if proper values are not entered in this data base.

Touching the blank field you want to fill out this window will appear:

A screenshot of a login dialog box. The dialog has a dark header bar with the word "Login" on the left and a close button (an 'X' icon) on the right. Below the header, there are two text input fields. The first is labeled "User:" and the second is labeled "Password:". At the bottom of the dialog, there are two buttons: "Cancel" and "OK".

Touch the blank field again an alphanumeric keyboard appears. Use it to write all the data. In this data base it is possible to set up the following values:

- Main motor in first speed: after this time the motor changes speed.
- Time ON pump first start: at the beginning of the cycle, a pump with deodorizer (when present) works for the time indicated.
- Time ON pump second start: During the cycle a pump with deodorizer (when present) starts again and works for the time indicated.
- Temperature delay OK: once that temperature is measured on the PLC, the temperature changes after the time indicated, in order to prevent that a noise signal has any effect on the cycle.
- Discharge door delayed opening: when the discharging temperature is reached PLC waits the time indicated before opening the hatch, to prevent vapour leaks.
- Product discharge time: the discharging phase will last the time indicated.

3.2.15. Setup sensor data


Fig. 45

In these submenus you can set up the sensor parameters. Not all values are possible, inside PLC some limits are configured and every value must respect these limits.

It is highly recommended to contact the producer before changing these values to avoid any malfunction during the cycle. Severe personal injury or substantial property damage can result if proper values are not entered in this data base.

Touching the blank field you want to fill out this window will appear:

Touch the blank field again an alphanumeric keyboard appears. Use it to write all the data.

3.2.16. Manual Function


Fig. 46

Fig. 47

In these submenus you can activate in manual some functions illustrated in par. 3.5

3.3. Switching on and switching off

Turn the "Main Switch" knob to position "I" or "ON".

Check for the presence of power supply (power light indicator on). If it is not, activate the differential knob of the machine.

At the beginning the CPU loads the internal program in 30 seconds, at the end the display shows this message:

EMERGENCY PRESSED OR MACHINE NOT STARTED

Verify that emergency stop buttons are not pressed, otherwise put them in the normal position rotating the red button, now push the button "START MACHINE" (the green light button must switch on), now the machine has started.

The display will show the same message above again, only after that the "START MACHINE" green light is on, press the icon ACK on the operator panel, the message will disappear.

If after this the machine will not start check the solving problem section.

In case of emergency the machine can be switched off by pushing one of the red emergency stop buttons.

Warning! Use this option only in emergency.

To switch off the machine press the "STOP MACHINE" button and turn the main switch on "0" or OFF.

Warning! During normal operation the "BY-PASS SAFETY DEVICES" light must be off, otherwise call authorized personnel in charge of the maintenance.

3.4. Automatic cycle

Put a bag under the discharging hatch in the steel collector and close it. With the cell full and the lid closed rotate the key selector to reach "Automatic" mode, wait 10 seconds until the machine closes the hatch (if the user forgot it open) and check the safety devices. Press the

button  . The automatic cycle will start.

If the machine does not start, check the display on the PLC panel for error messages by pressing F2. Check the display and follow the instruction in paragraphs 3.6 and 3.7.

Cycle can begin only if the display shows

NO MESSAGE

The signal for the end of the operating cycle is given by a red light located above the main control panel, which gives both a light and sound signal. To halt this signal, turn the mode key to reach "Manual" position. When the cycle is finished and the bagging machine has completely stopped (where this accessory is present, before taking out the bag it is important to shake the filter of the bagging machine so all the material attached to the filter can drop inside the bag).

Now the machine is ready for another cycle.


WARNING! Always wear anti-dust mask to close and change the sterilized material bags.

3.5. Manual machine controls

To use the 'Manual' functioning mode, all you have to do is rotate the selection mode key on "MANUAL" position.

By pressing buttons on the keyboard the following functions can be carried out:

| | |
|---|--|
|  | Start the automatic cycle |
|  | LED on - start motor at low speed, LED off - STOP motor |
|  | LED on - cabinet ventilator filters in operation, LED off - ventilator stopped |
|  | LED on - resistance operating, LED off - resistance off |
|  | LED on - general water valve open, LED off - valve closed |
|  | LED on - water valve in cell open, LED off - valve closed |
|  | LED on - door in opening phase, LED off - door closed |
|  | LED on - door in closing phase, LED off - door closed |
|  | LED on - lid unlocked, LED off - lid locked |
|  | Increases the speed motor |
|  | Decreases the speed motor |
|  | LED on - motor turns in clockwise, LED off - motor turns in counterclockwise |
|  | Switch on the main motor servo fan |
|  | LED on - pump is working, LED off - pump stops |

Note that when the LED is green it means that the function is on in manual mode, while when the LED is blue it this means that the function is on in automatic mode. The normal operation of the machine will be in 'Automatic' mode.

- The operations carried out in 'Manual' mode are exclusively for
- checking the functioning of the system or of other components when empty or when loaded,
 - solving problems which occur during the execution of the automatic cycle,
 - carrying out emptying operations and cleaning the sterilization cell.


3.6. Service messages

To access the level of alarm signalling, press the button F2 or the icon. In this way you enter the alarms menu and by scrolling through the messages using the arrows you will find that any alarms present are shown on the display.

If the display shows this message.

HANDLE OPEN


It means that the handle is open, in this case the automatic cycle cannot start and in manual

mode the motor cannot rotate. Push the key  (verify the key has the led on), close the handle and press the button again, this time the led must be off, force the handle a little. If the message is still showing, contact authorized technical assistance.

If the display shows this message.

LID OPEN

In this situation the automatic cycle cannot start and in manual mode the motor cannot rotate

because the cover is not closed. Even in this case push the key  and verify that the led is on, take the handle of the cover and close it. The message should disappear.

If the message is still on contact authorized technical assistance.

If the display shows this message:

DISCHARGING HATCH NOT BLOCKED

The discharging hatch is not closed correctly, the automatic cycle cannot start. Put the

machine in manual mode, open the discharging hatch pushing  clean the surface of

the hatch well and close again with the key . The message should disappear, otherwise contact the technical assistance.

When pushing ESC the display shows this message.

NO MESSAGES

Now the cycle can start normally.

3.7. Alarms

In this paragraph the main alarm messages and the procedure for resetting the machine for normal operation after an alarm message are examined.

To access the level of alarm signalling, press the F2 button. In this way you enter the alarms menu and by scrolling through the pages using the arrows on the screen you will find that any alarms present are shown on the display.

If during the normal operation of an automatic cycle the PLC detects a reason for an alarm, the cycle will immediately be interrupted and the display will begin to flash indicating the alarm, a red light will go on in the alarm column and you will also hear an acoustic signal. To reset the correct functioning of the machine, follow one of the next procedures, according to the nature of the alarm

LIST OF ALARMS

Immediately after turning on the general power switch and after having pressed EMERGENCY STOP buttons, the screen displays the following message.

EMERGENCY PRESSED OR MACHINE NOT STARTED

Verify that the two "**EMERGENCY STOP**" buttons are not pressed, otherwise replace them to the normal position by triggering them off with counterclockwise turn.

Press "START MACHINE", the button should light up, now press the ACK icon and the message will disappear.

If you press the START MACHINE button and the light does not turn on verify that the "**EMERGENCY STOP**" button is not pressed and check that the protection panels of the machine are tightly screwed in. Press START MACHINE, if the machine still does not start then call technical assistance.

If the display shows

MAXIMUM CURRENT MAIN MOTOR PROTECTION

turn off the power, positioning the main switch in the "0" (off) position, open the inverter group board and find the breaker marked MT5 and push the blue button to reset him.


Fig. 48

Then close the inverter group board and turn the main switch to the "1" (on) position. Finally, follow the instructions at the bottom of this list of alarms.

If the display shows

THERMAL PROTECTION FUMES SUCTION

turn off the power, positioning the main switch in the "0" (off) position, open the electric board and find the breaker marked MT2 and move its lever to the "1" (on) position. Then close the switchboard and turn the main switch in the "1" (on) position. Finally, follow the instructions at the bottom of this list of alarms.

If the display shows

THERMAL PROTECTION FILTER SUCTION MOTOR

turn off the power, positioning the main switch in the "0" (off) position, open the electric board and find the breaker marked MT1 and move its lever to the "1" (on) position. Then close the electric board and turn the main switch to the "1" (on) position. Finally, follow the instructions at the bottom of this list of alarms.

If the display shows

CHECK MT11 PROTECTION OF MAINMOTOR SERVO FAN

turn off the power, positioning the main switch in the "0" (off) position, open the electric board and find the breaker marked MT11 and move its lever to the "1" (on) position. Then close the electric board and turn the main switch to the "1" (on) position. Finally, follow the instructions at the bottom of this list of alarms.

If the display shows

LOW WATER PRESSURE

check that both the water tap is open and that there are no obstructions along the hose that carries the water to the machine. Then follow the instructions at the bottom of this list of alarms.

If the display shows

WATER PRESSURE RECIRCULATING SYSTEM TOO LOW

please check that the water recirculating system pump is working.

If the water recirculating system pump is working first of all you have to check if there is enough water in the circuit and, if necessary, refill with other water.

If the problem remains contact the assistance service.

If the display shows

CHECK OIL FILTER AND CLEAN IT

it is necessary to clean the filter, shown in red in the picture, of the oil from the machine.


Fig. 49

To do it it is necessary to unscrew the N°6 M6 and open oil pump vessel. Then follow the instructions at the bottom of this list of alarms.

If the display shows

OIL PRESSURE TOO LOW

please check that the oil pump is working; if not probably the breaker MT7 switch off. In this case turn off the power, positioning the main switch in the "0" (off) position, open the electric board and find the breaker marked MT7 and move its lever to the "1" (on) position. Then close the electric board and turn the main switch to the "1" (on) position. Finally, follow the instructions at the bottom of this list of alarms.

If the oil pump is working first of all you have to check if there is enough oil in the circuit. To do it unscrew the cap, shown in blue in the picture, and, if necessary, refill with oil of the same type.

If the oil pump is working and if there is enough oil please clean the the filter, shown in green in the picture, of the oil to the machine.

If the problem remains contact the assistance service.

If the display shows

CHECK MT9 AND MT8 THERMAL PROTECTION DISCHARGING PUMP AND COOLER

turn off the power, positioning the main switch in the "0" (off) position, open the electric board and find the breaker marked MT9 and MT8 and move its lever to the "1" (on) position. Then close the electric board and turn the main switch to the "1" (on) position. Finally, follow the instructions at the bottom of this list of alarms.

If the display shows

OVERLOAD INVERTER

it means that a fault stops the drive. You must reset the drive and find a solution to the problem.

Before you contact the distributor or the factory because of unusual operation, prepare some data. Write down all the texts on the display, the fault code, the fault ID, the source info, the Active Faults list and the Fault History.

To reset the drive switch off the main switch for 60 seconds.

If display shows

TEMPERATURE SENSOR BROKEN

it means that in 22 minutes the temperature in the cell do not rise up to 40 °C. If you are at the first cycle of the day and that the machine was cold it is possible that the sensor is not damaged and it is only necessary to restart the cycle.

If the problem remains contact the assistance service.

If the display shows

MOTOR'S WINDING TEMPERATURE TOO HIGH

it is necessary to finish the cycle and to empty the cell. Now start the machine in manual mode and let the motor works without material until the motor temperature display shows less than 85 °C.

If display shows

CELL TEMPERATURE TOO HIGH DO NOT OPEN THE LID

it means that the cell's internal temperature is excessively high, at a level which could damage its external covering. In this case you need to wait until the temperature shown on the indicator falls to 90°C before switching on the machine again (follow the instructions at the bottom of this list of alarms).

WARNING. For no reason open the lid and pour in water in manually.

Warning! This alarm could be due to a malfunction in the temperature probe. Before using the machine again in the normal way we advise that you consult the appropriate paragraph in the chapter regarding maintenance and, if necessary, contact the assistance service.

If display shows:

INSPECTION DOOR OPEN

check the inspection door and if it is necessary, close them. Then follow the instructions at the bottom of this list of alarms.

Once the cause of the alarm has been eliminated, the display continues to flash until the


alarm has been acknowledged. To acknowledge the alarm, press icon . If the LED is lit because an alarm is indicated, pressing this key will switch it off, and cause the panel cursor to stop flashing.


However, if other alarms are also present, the display will continue to flash to indicate them, in which case repeat the sequence just described until all have been removed.

If after this procedure the display shows the alarm again, please contact authorized technical assistance.

3.8. Filling of sterilization vessel

Put the machine in manual mode, open the hatch by pressing  and clean the metallic surface around the discharge door with the special device supplied with the machine, then


close the hatch again by pressing the button , now press the button  (a led button lights up). With this operation the lid handle will be released. Open the lid and then


press the button  again (this time the led light must go off) so the lid is blocked, now filling the vessel with the potentially infectious waste.

When filling the vessel, the lighter bags must be placed above the blades to facilitate the pulverization process, while the heavier ones must be placed last. If possible fill the cell with material from different wards.

The vessel must be completely filled but without exerting pressure on waste, do not compress waste.

Please remember that the machine is for solid waste and the cycle cannot finish with humidity higher than 40%.

After the vessel has been filled with wastes, press  (the led must be on) and close the

lid and handle, now press again  (the LED must be off) and force the handle forward. Now the machine is ready for a cycle.

Do not fill the vessel with the discharging hatch open.

Warning! The lid surface can be very warm. Grip the lid only using the plastic handles.

Do not fill the vessel with lid unblocked. The LED on the key  must be off.

3.9. Carry out a pseudo-cycle in 'Manual'

Although, as previously mentioned in the present manual, the normal functioning of the machine is in 'Automatic' mode, it is also possible to carry out a pseudo-cycle in 'Manual' mode to check the machine or finish a failed cycle.

To do this, first of all load the material to be treated. Make sure that the loading door is closed, and secure the lid tightly.

Put the machine in manual mode.

Start the rotor by pressing the  key, and at the same time start up the fan of filters by pressing , the resistance  and the pump .

After 60 seconds, press the  key once again. Now the green LED will be off and the pump will stop working.

When the temperature reaches 50 °C, press .

Observe the ammeter. Once the consumption of electric current has fallen below 100 A, wait 2 minutes while checking that the ammeter does not exceed this threshold. If it does not,


increase the speed motor by pressing the  key until the inverter gains 5 Hz. If the ammeter does not exceed 100 A repeat this operation after 2 minutes.

Observe the ammeter. If the ammeter does exceed 120 A decrease the speed motor by pressing the  key until the inverter loses 5 Hz.

Leave the rotor to rotate until the temperature of 150°C is reached. At this point press the  key until in the operator panel you can read 20 Hz and the  key in order to switch off the resistance.

After having reduced the motor speed, press the  button and cool the mass until the temperature falls to about 95 °C. Press the  button again in order to stop the cooling and wait a minute in order to allow the water to evaporate. After stopping the "water in cell", please check if "general water" is activated, otherwise activate it by pressing .

At this point it is possible to unload the treated material by pressing the  button. The unloading of the material will take about 2 minutes.

Once the unloading has finished, stop the rotor by pressing the  button.

Allow 2 minutes to pass, then close the unloading hatch by pressing the  key.

Now press the **PUMP** button and operate the pump for 30 seconds.

This operation causes disinfectant to be sprayed into the cell.

Wait 1 minute to let the spray be completely sucked into the absorption column, then press

and

to switch off the ventilator and the electrical valve that pours water in cell.

From this point on it is possible to open the lid and start loading operations for a new cycle.

3.10. Downloading material from the cell in MANUAL mode.

Sometimes the automatic cycle can be interrupted after reaching a temperature of 150° C, in this case it is not necessary to begin the cycle again as the material is already sterile, it is sufficient to empty the vessel using the following procedure.

Make sure the internal temperature of the vessel is inferior to 95° C, that a bin with an empty bag is placed under discharge hatch and the lid for the discharge inspection is in position.

Start the motor by pressing , set the speed at 20 Hz pressing the  button and

after 10 seconds open the discharge hatch , wait 2 minutes then press the  button, the motor should turn off. Now the cell is empty you can start a new cycle.

Never discharge material with the motor on speed higher than 20 Hz.

Never open the lid to inspect material when the motor is on, fragments of sterile waste may be violently thrown out.

3.11. Verifying the Effectiveness

The sterilization cycle of the machine can be verified by using biological indicators.

Inside the vessel (Fig. 5) there are three compartment for the biological indicator, before filling the vessel unscrew the cap (A in Fig. 5) with a 17 mm hexagonal key and insert a testing vial. Be sure that the closing cap is properly tight to avoid the possible breakage of the vial.

Always use leather gloves as the inside of the vessel remains hot even after the cycle. Refill the vessel and perform an automatic cycle. At the end of the automatic cycle open the lid and using leather gloves and an hexagonal key unscrew the cap again and extract the vial with caution. The vial is available for sterilization verification.

Fig. 50

WARNING! Only use glass bioindicators with diameters not more than 11 mm.

WARNING! The NEWSTER®NW30 sterilizer does not work with high pressure, therefore the effectiveness cannot be verified by using biological indicators that are activated only by pressures above atmospheric pressure.

3.12. Integrated Waste Collection system

The bagging system is integrated in NEWSTER®NW30. The sterilized material is automatically discharged into a steel box connected by an airtight seal to the discharge system.

At the end of the cycle, after the discharging phase is completed, the steel box is take away using the support with wheels.

In case of trouble with the discharging system an alarm appears on the control panel (see paragraph 3.6).

It is necessary to replace the sack for every cycle.

Fig. 51 - The automatic bagging system

The discharging sequence:

Fig. 52

Fig. 53

Fig. 54

4. Maintenance

4.1. Informative notes

The complete maintenance of NEWSTER®NW30 machines cannot be covered in these paragraphs. A thorough maintenance can only be carried out by qualified personnel who have attended a training course by the producers or the authorized distributors.

4.2. General norms

Maintenance must be carried out by specialized technicians who have been trained in the following specific areas:

- Mechanical maintenance
- Electrical maintenance

The professionalism and competence of the technician, lies with the person responsible for safety.

Before starting maintenance intervention the person in charge of safety must:

- Vacate the work area of those not involved.
- Make sure that the necessary instruments are at hand and in good condition.
- Verify that there is sufficient light, otherwise provide portable 24 volt lights.
- Make sure the authorized personnel are equipped with the safety regulation clothes (protective gloves, goggles etc).
- Verify that the maintenance technician has carefully read the instruction manual and understands the function of the machine.

Before starting the technician must disconnect all supplies (electrical and pneumatic) and ensure the security block on the machine is in function.

Warning! If it is absolutely necessary to work on the machine during the maintenance, the service engineer must keep at a safe distance and the safety blocking switch must be close at hand.

After the intervention and before restarting the service engineer must check the security system and the complete performance.

Work on the electrical motors must be executed ONLY by the service engineers who have been trained by either the producer or the distributor of the machine.

Warning! After every service the engineer is responsible for the complete check and safety of the machine and the protective operative system

Maximum dependability of the machine at a minimal cost is the result of a maintenance program and a careful inspection during the life of the machine

Warning! Before starting any maintenance or cleaning service the machines energy sources must be shut down, meaning the disconnecting switch and the air compressors loading valve.

Warning! Before starting any maintenance or cleaning service the engineer must be wearing the protective clothing, goggles and gloves.

Warning! Before starting any type of service on the machine affix a sign on the machine indicating:

**MACHINE MAINTENANCE
DO NOT ADD MATERIAL**

Warning! The mechanical instruments inside the cell are very sharp, therefore always use protective gloves.

Before the re-starting operation of the machine, check entire system in accordance to the start up program. After every service carry out a few trial cycles thus verifying the correct function of the machine.

If these safety precautions are not followed severe harm may be caused to personnel.

Before any work:

- Check that the correct label is on the machine.
- Check that the water tap of the machine is open.
- Check the water discharge tube of the machine.
- Check if any tools were left inside the machine.

Any intervention to the sterilizer by person or persons not qualified or not authorized by producer or distributor void the guarantee.

Use original spare parts or parts approved by producer or distributor.

4.3. Replacing consumable material

This paragraph will examine the material that the operator needs to replace. For maintenance interventions that are not treated in this manual please contact your after sales service.


Any intervention made on the sterilizer by unqualified or unauthorised personnel will make any warranty made on the sale of the machine void.

4.3.1. List and duration of consumable parts

In the following tables you will find the stimated duration of consumable parts.

Due to a stronger use of the machine it is possible a shorter duration of the parts. Please check, at the end of each work session, the consumption of the consumable parts.

GRUPPO LAME-BLADES GROUP, GRNW050-03-03

Fig. 55

| BLADES GROUP, GRNW050-03-03 | |
|------------------------------------|------------------------------------|
| CODE | DESCRIPTION |
| FE00195 | Self-Locking nut M10 zinc metal |
| FE00226 | Screw TCEI Zinc M10x25 |
| FX00141 | Screw TCEI Inox M8x50 |
| FX00226 | Hex Screw inox M20x65 |
| FX00260 | Screw TCEI Inox M8x20 |
| FX00263 | Screw TCEI Inox M12x30 |
| FX00276 | Screw TCEI Inox M10x70 |
| FX00283 | Screw TCEI Inox M6x16 |
| FX00287 | Washer inox 6,5x18x1,5 |
| FX00296 | Lock washer tooth for M6 |
| M000829 | 2 Holes Blade |
| M100830 | Head Biade |
| ME01551 | Support vertical blades |
| MC01553 | Sheet srew Lock, key 30 |
| 11101557A | Fixed fLange Labyrinth |
| M101558 | Joint Plate Blades Support Holder |
| M101559 | Moving flange Labyrinth |
| M101614 | Radial shaft seal, 85x110x12 Viton |
| ME01616 | OR 4625 158,3 X 3,53 |
| M101617 | OR 6775 196,21X5,34 |
| MC01626 | Tongue, section 14x9 L=40 |
| M101761 | Tongue, section 14x9 L=25 |
| M101812 | Blades support, main part milled |
| MF01552A | Single Vertical Blade NW50 |


Fig. 56

| CELL GROUP, GRNW050-04-02 | |
|----------------------------------|---------------------------------------|
| CODE | DESCRIPTION |
| A500279 | Piston 0.20 C.15 Single Effect |
| FC00199 | Washer Zn 10.5x30x5 |
| FX00229 | Washer 24,5x39x3 Stainless |
| FX00261 | Screw TCEI Inox M8x25 |
| FX00450 | Nut M15 Short Stainless |
| FX00452 | Nut M24 Short Stainless |
| FX00455 | Washer 17x23x3 Stainless |
| FX00457 | M10 NUT A2 SHORT |
| M000289 | Gas Spring Mod. 198.390.35 |
| M000391 | Seal Edge Cell |
| ME00456 | P1100 Probe |
| M100517 | Vassel Sleeve Sp.35 |
| M000739 | Cover Short Sensor 35 |
| MC00751 | Cover Long Sensor Vassel 35 |
| MC00752 | Protection measuring hole Vassel 35mm |
| MC01588A | Protection Vassel NW50 |
| MC01588B | Protection Vassel NW50, Sensor Holden |
| MC01612 | V-Ring Seal V-110A |
| MC01658 | Vessel Protection 2nd lever, Hardox |
| MC01659 | Discharge Mouth Protection NW50 |
| MC01663A | Bottom Vassel Protection Hardox NW50 |
| MF00742 | Sensor Short Temperature |
| SANW050-02-02 | Lid NW50 |


Fig. 57

| SUBASSEMBLY NW50, SANW050-02-02 | |
|--|-----------------------------------|
| CODE | DESCRIPTION |
| A000026 | Coupling 6 mm for Air Hose m 1/4' |
| A000031 | Bend M/F 1/4' 90° |
| AC00124 | Bulkhead G1/4' M20x1,5 |
| A000173 | Atomizer Nozzle |
| A000505 | Hose 6x4 Neutral PTFE (15cm) |
| AC00505 | Hose 6x4 NeutrIL PTFE (15cm) |
| FX00230 | Washer Inox 8,5x32x2 |
| FX00260 | Screw TCEI Inox M8x20 |
| FX00280 | Washer 13x24x2,5 Stainless |
| MC00676 | INSPECTION VESSEL LID |
| MC01611 | LID GASKET L=2560 |
| MC01670 | REAR FRAME LID |
| MC01671 | LID FRAME |
| ME05020 | Nozzle Clamp |
| ME05024 | Nozzle Sheet |
| MF01666 | CIRCULAR LID |


Fig. 58

| GRNW 010-09-03 | |
|-----------------------|----------------------------------|
| CODE | DESCRIPTION |
| MC00369 | Absolut Filter for Newster NW 10 |
| MC00368 | Filter, Active Carbon |


Fig. 59

| SANW 050-01-00 | |
|----------------|--------------------------|
| CODE | DESCRIPTION |
| MC00157A | Net, Filled Bodies Block |
| MC00156 | Filled Bodies |
| MC0180 | INOX Filter |
| MC01578A1 | Net, Filled Bodies Block |



Fig. 60

| GRNW050-21-01 | |
|---------------|-------------------------------|
| CODE | DESCRIPTION |
| MF00841 | Heat exchanger NW50 assembled |
| AC00093 | Pipe 12x10 blue Rilsan |
| MC00456 | PT100 Probe |

S.A. ACQUA-WATER SUBASSEMBLY, SANW050-24-00

Fig. 61

| SANW050-24-00 | |
|---------------|------------------------|
| CODE | DESCRIPTION |
| MF00843 | Heat Exchanger NW50 |
| MF00842 | Exchanger below part |
| FX00287 | Washer Inox 6,5x18x1,5 |
| MF00844 | Exchanger top part |
| FX00285 | Nut M6 Stainless |
| FX00220 | Grower Washer Zn M6 |
| MC00547 | Spacer Bush |
| FX00295 | Screw TCEI Inox M6x20 |
| MC01741 | OR CORDA 250,2 x 5,34 |



Fig. 62

| GRNW 050-06-01 | |
|----------------|----------------------------|
| CODE | DESCRIPTION |
| MC01548 | Plate Blocking Gasket Door |
| MC01535 | Seal Discharging Hatch |

4.3.2. Changing printer paper

Use only thermal printer paper. Width 57,mm (+/-0,5mm) and the diameter should be less than 50mm. To change paper do as follows (see diagram below):

- Open the printer's door and press on the support of the tilting mechanism of the printer where it indicates PUSH:
- Place the roll of paper, respecting the direction of rotation as shown in figure:
- Press the FEED button to let a few centimeters of paper out of the printer and push again on the support mechanism of the printer where it indicates PUSH, then replace in the original position:
- Rip the paper and close the door by inserting the card slot in the door itself.

4.3.3. Vertical blade changing.

Vertical blades (MF01552) cut and open wastes bags, to spread the wastes inside the cell. They should be changed, when problems with opening bags appear or after 1000 cycles. For the substitution of the blades please follow the described procedure:


Fig. 63

- With a N°6 Allen key unscrew the FX00283 bolt
- Remove the MC01553 component.

- With a 24mm wrench mounted on an extension unscrew bolt FX00226.
- With Allen key and an extension take off the two FX00268 screws.
- Take blade support from cell. Do it wearing protective gloves to avoid injuries.
- Also take off MC01551.
- Clean MC01551 and all contact surfaces.
- Check all screws. Change if it's necessary.
- Install new vertical blades. Use only original parts.
- Place the parts on there place, tight the screws.
-

Warning. Parts installed inside cell have sharp edges. Always use protective gloves when working.

4.3.4. Blade support dismanteling.

Blade support (MC01550) handle with strong dynamic events. Specially for this, it is made with special hard and resistant steel. It should be changed every 1500 cycles..


Fig. 64

- Dismount vertical blades, as is described in chapter 4.3.3.
- Unscrew with hexagonal key on extension bolts 221.
- Rise the blade support from cell. Use protecting gloves.
- Clean well all contact surfaces before installation of new blade support.
- Mount all parts in their places. Tighten the screw with the correct torque.

Use only original or authorized spare parts. Using non-original or non-authorized parts may cause malfunction of the machine, or accident.

Warning. Parts installed inside cell have sharp edges. Always use protective gloves when working.

4.3.5. Changing blades.

Fragmentation and heating of the material is made by the blades. For that reason, to keep maximum efficiency of the machine, it is important to pay extra attention to the condition of the blades. Change them every 200 cycles.

- Dismount the blade support from the machine, and place it in strong vise.
- Use an N°8 Allen key and a N°17 wrench to unscrew all 12 bolts. The bolts are locked by nuts placed on the bottom side of blade support.
- Clean contact surfaces of the blade support and blades.
- Change used bolts with new ones.
- Place new blades and tighten the screw.

Fig. 65

Warning. Parts installed inside the cell have sharp edges. Always use protecting gloves when You work.

4.3.6. Replacing the short temperature sensor.

The short temperature sensor (called hereafter short sensor) is responsible for the temperature monitoring signal which is transmitted to the PLC and printer, to control the cycle.

The machine is designed in a way to protect short sensor. But it can happen, that replacing will be needed. To do it, follow next procedure.

- Switch off the machine and take off the panels.
- Disconnect the two wires connected to the short sensor. (Fig. 64).
- With a N° 7 and a N° 17 wrench, release locking (Fig. 65), it is not needed to unscrew it completely, it is enough to loosen it.
- With a N° 7 wrench unscrew the short sensor (Fig. 66).
- Place a new sensor, make sure that it is completely screwed in place.
- Lock the sized 10 nut with a N° 17 wrench holding the short sensor with a N° 7 wrench.
- Connect wires to the sensor. The order of connection is not important.

- Close the panel of the machine. Make sure that protection circuit is active.

Fig. 66

Fig. 67

Fig. 68

The short sensor is one of the parts responsible for the correct work of the machine. Use only sensors supplied by the producer of the machine, or authorized sellers.

4.3.7. Replacing the short sensor cover.

As explained in paragraph 4.3.5 it is very important to protect the short sensor, so please pay special attention to the condition of the short sensor cover. Replace the sensor cover every 250 cycles.

- Remove the short sensor as described in the previous paragraph, to reach the situation shown in Fig. 66.
- With a N° 24 wrench remove the locking nut holding the main nut (Fig. 67).
- Unscrew the main nut with a N° 24 wrench, holding the sensor cover with a N° 13 wrench (Fig. 68).
- With a rubber hammer hit the sensor cover, to put it inside cell (Fig. 72).
- Clean the contact surfaces of the sensor cover and cell well.
- Use red silicon (temperature resistance till 250°C) to seal the connection between the sensor cover and the cell.
- Replace the sensor cover with a new one, mount all parts in place, be sure that all the nuts are tight.
-

Fig. 69

Fig. 70

Fig. 71

Fig. 72

WARNING! The sensor cover has influence on the temperature measurement. For this reason use only original spare parts.

4.3.8. Replacing the long temperature sensor.

The long temperature sensor is part of the protection circuit. For this reason use only original spare parts, to avoid malfunction or risky situations during use of the machine.

- Release the plastic cover and pull out sensor.
- Disconnect wires connected to the temperature display.
- Install a new sensor. The polarity is important, so take care with the connection. Watch carefully how the previous sensor was connected, and connect the new one in that same way.
- Tighten the plastic cover.

4.3.9. Replacing the long temperature sensor cover.

Every 250 cycles replace the cover of the long temperature sensor, following the procedure described below.

- Remove the long temperature sensor, as described in the previous paragraph.

- Unscrew the sized 10 nut with a N° 17 wrench holding the sensor cover with a N°8 wrench (fig. 73).
- With a rubber hammer hit the sensor cover to insert it inside the cell (fig.74).
- Clean the contact surface of the cover.
- Place a new cover, and put all parts in places. Tighten all the nuts well.

Fig. 73

Fig. 74

The long temperature sensor is a part of the protection circuit. For this reason use only original spare parts, to avoid malfunction or risky situations during use of machine.

4.3.10. Changing the lid gasket.

Every 600 cycles it is advised to change the gasket of the lid (MC01661), following this procedure:

- Remove the old gasket from it's seat.
- Insert a new gasket beginning from point A (Fig. 75).
- To place the new gasket in position you can use a screwdriver. Be careful not to puncture or cut the gasket.
- When the gasket is in place leave an extra 3 cm of lenght and insert one end into the other.


Fig. 75

4.3.11. Replacing the carbon and absolute filters.

Carbon filters, which are placed inside the filter group should be changed every 6 months. The absolute filter should be changed every year.

To replace the carbon filter follow the procedure described below:


Fig. 76

To replace absolute filter follow the procedure described below:


Fig. 77

WARNING! Filters are part of the air treatment circuit of the machine. Use only original or authorized parts.

4.3.12. Maintenance of deodorant system

The deodorant system located on the lid consists of an atomizer fed by compressed air and enzymatic deodorant liquid.


Atomizer

Fig. 78

Every 200 cycles proceed with the atomizer cleaning in order to prevent any pipes clogging. Proceed with the following steps:

1. Remove the three screws and the MC05024 plate.

2. Unscrew the nut with a hex key 22 and remove the AC00173C e AC00173D nozzles, taking care not to lose the AC00173B seal.

3. Clean up the AC00173C and AC00173D nozzles, possibly with compressed air.
4. Reassemble the components.

Daily check the level of the enzymatic deodorant in the tank inside the stainless steel container placed next to the main electrical panel.

The consumption is about 0,25 liters per cycle. Therefore, every 12 cycles it is necessary to fill the tank with a solution composed of 1/8 enzymatic deodorant (Activair) and 7/8 tap water.

4.3.1. Checking and replacement of the cartridge filter

A cartridge filter, to protect the system from any impurities contained in water, is placed at the pump outlet.


Fig. 79

Every six months it is recommended to open the filter by unscrewing the blue cap and visually check the status of the cartridge. Wash it if dirty or replace it if damaged.

4.3.2. Add new glycol

Fig. 80

To avoid water freezing problem, it is recommended to annually fill the water tank with a proper amount of ethylene glycol or propylene glycol solution (an example of glycol product is shown in figure 39 – Newster code: MC00130). Open the tank by unscrewing the black cap at the top of it and pour the solution inside it. To add the right amount of solution, follow the instructions on the product data sheet. Close the cap.

| Ethylene glycol | |
|---------------------|---------------|
| Freezing point [°C] | |
| X [vol %] | T freeze [°C] |
| 10 | -3.6 |
| 15 | -6.1 |
| 20 | -8.9 |
| 25 | -12.1 |
| 30 | -15.7 |
| 35 | -19.9 |
| 40 | -24.8 |
| 45 | -30.4 |
| 50 | -37.1 |
| 55 | -44.5 |
| 60 | -50.5 |

| Propylene glycol | |
|---------------------|---------------|
| Freezing point [°C] | |
| X [vol %] | T freeze [°C] |
| 10 | -3.4 |
| 15 | -5.3 |
| 20 | -7.4 |
| 25 | -10.0 |
| 30 | -13.1 |
| 35 | -16.9 |
| 40 | -21.5 |
| 45 | -27.0 |
| 50 | -33.7 |
| 55 | -41.6 |
| 60 | -50.9 |

4.4. At the end of the cycle

At the end of every cycle you must clean the machine.

Using a thoroughly clean wire brush clean:

- the mesh under the lid (close the port that leads to the condensation column);
- the edges of the discharge port and hatch.

The machine will now be ready to start another cycle.

End of the work session

Everyday at the end of the shift the operator must carry out simple operations that serve to keep the machine clean and efficient. Following are the operations that must be carried out at the end of each shift when the machine is turned off:

- Clean the metal mesht situated on the lid with a vacuum cleaner, if necessary rub the net with a brush and remove the dust on it:


Fig. 81

- Clean the bottom of the cell controlling wear of the blades:
- Check the inside of the cell making sure the screws are tight:
- Clean dust from the nozzles that are located within the cell:
- Disconnect the tube that leads vapours from the machine group filters and empty if condensation has formed during operation.
- Wash the sterilization vessel. Repeat this operation twice

Leaving the machine unattended

At the end of the working day it is advisable to leave the machine with the lid and the handle locked in place. After having closed them, it is necessary to turn the machine off, move the main switch to the "0" position and use the key to lock the door of the control panel. The key must be kept by the operator.

4.5. Maintenance for long inactivity

If you are planning on not using the machine for long periods of time procede as follows:

- Switch off the general power on the electrical board.

- Thoroughly clean the cell and oil with rust proof oil the bottom of the cell and the rotor blades.
- Clean the filter located on the lid with a metal brush.
- Immerse the hypochlorite sodium filter in water and wash the tubes rinsing well with water.
- After washing remove the filter from the water and leave it in open air, activate the pump until the circuit is completely empty.
- Empty loading and unloading water tubes.
- Empty water present in the tube that connects the machine to the filter cabinet.
- Protect the machine and all its components from dust.

Finally, in order to keep the machine in perfect working efficiency, even during long periods of non use, we advise that you contact specialised personnel.

5. Trouble shooting

5.1. Rotor blocked because of overload

If during the cycle the rotor stops because of an overload, proceed as described in paragraph

3.5 . Now push the  key (The LED will light up and the motor will turn in clockwise) and restart the cycle.

Remember that during the normal functioning it is preferable to maintain the motor in counterclockwise direction (LED on the  key OFF). Use position clockwise direction only to superate conditions blocking the rotar .

There is the possibility with certain malfunctions that the cycle cannot re-start .

If it is not possible to re-start and complete the sterilization cycle , the material and the equipment can present biological risks .

Place the machine in "Manual" mode. Open the cell's lid and the discharge hatch, empty the machine into a few bags by pushing the waste through hatch. Load the bags in a new cycle to be reworked.

In the event that a malfunction occurs right after the start of the cycle and the waste has not been sufficiently shredded to be able to extract from the bottom hatch, proceed to extract it by hand from the top after protecting oneself with a smock, heavy-duty gloves and a mask against particulate matter (follow the hygiene regulations closely). At the end of the operation the smock must be washed immediately in water and disinfectant and the other safety equipment must be bagged and treated like hospital waste.

Remember that , by loading big bags containing material that is particularly resistant like, for example sheets, scrubs, rubber and plastic tubes, and if at the same time the blades are not very sharp, it is possible that the resistant torque vastly exceeds that available from the motor. In that event , in fact, the rotor instead of penetrating the material drags the bags forward without shredding them resulting in overload.

5.2. The printer does not print and the cycle does non finish

If during the cycle : the printer does not print any values , the temperature shown on the display of the operator's panel doesn't rise and at the same time the temperature display on the machine's console shows that the temperature is rising ; the temperature probe is broken. Replace both the temperature sensor and the sensor cover as described in the maintenance chapter.

If the printer still does not print the temperature values after the replacement the cause may be the interface between the sensor and the PLC . In this case call specialized assistance.

5.3. The treated material is too hot when discharged

The temperature allowed for the discharge of the material is 95°C.

If the material is discharged at a temperature of over 120°C control the system introducing water.

Especially the water faucet , the injection nozzles inside the cell and the solenoid valves EVO and EV1 .

Check the state of wear of the short temperature sensor which commands the PLC.

If after these controls the material is discharged at a temperature above 100°C contact technical assistance.

Attention ! Keep in mind that material that is too hot involves high risks of spontaneous combustion even several hours after being discharged and bagged.

5.4. Strange noise of metal banging inside the cell

The machine is not designed to treat large pieces of metal.

If there are hard, compact metal parts with a weight of over 50 grams, the machine will produce a strange noise caused by the banging of metal pieces against the rotor and the sides of the cell.

In these cases, to reduce as much as possible, the possibility of damaging the rotating and fixed blades , it is necessary to slow down the speed of the rotor ending the cycle in "Manual" mode in low speed.

To do so it is necessary to shut down and restart the cycle in "Manual" mode at low speeds. Complete the treatment by performing all operations (heating, cooling, exhaust) as described in section 3.6, but do not go into high speed. Of course, the duration of the cycle will be longer . At the end of processing open the cell and ensure that there are no more metal masses on the bottom of the cell.

5.5. The machine heats up less than usual

Note that if the cell has been filled and the rotor turns, it is impossible that the machine does not reduce heat.

If the temperature remains practically constant and it seems that the warming is weak , check the status of the temperature sensor and of the blades; if the blades are very worn their warming effect is greatly reduced .

It is also possible that the blades are a little worn but the material is packed in bags so light that they "float" and do not descend into contact with the rotor. In this case the material in contact with the rotor is very hot but the average temperature inside the cell remains fairly low. Normally, in these cases , before the end of the cycle the material above reaches the rotor and the temperature rises very quickly. If this phenomenon occurs frequently, the lighter packages must be loaded into the cell first, followed by the heavier ones

However, if these were not the cause , then the fact is due to:

- load of waste is too light
- moisture content of the waste is too high
- Remember that the machine is a solid waste sterilizer , so the weight of the material must not be less than 10kg and the percentage of humidity should not exceed 40-45%

5.6. Steam escapes from the cell

The output conduit of the vapors from the cell must remain free so that the suction due to the fan of the filters group maintains the cell in slight depression. In this way, the vapors can not get out of the cell by any other way except through the filter group.

If during the cycle the depression fails, a leak of vapours from the seal of the cell may occur. It is possible that excessive dust has clogged the circular filter on the lid; in this case it is necessary to open the cell and thoroughly clean the filter and then check that the vapor suction hole is not blocked.


Fig. 82

If even after this intervention the steam continues to exit from the cell, put the machine in "Manual" mode and turn on the vacuum cleaner filters in the group, cut a strip of paper 4 cm wide and 20 long and bring it close to the vapour suction hole holding it tightly ; if the paper is not sucked into the tube it is sure that the columns of lowering temperatures are too clogged with limestone. In this case the intervention of specialized technicians is required , please contact support.

5.7. Table of the main problems, their causes and solutions.

| PROBLEM | PROBABLE CAUSES | SOLUTION |
|---|--|--|
| Loud squeal coming from the transmission group | Belt not stretched | Tighten the belts and check their condition, change the whole group, if necessary (MAINTENANCE SPECIALIST) |
| Abnormal vibrations of the machine running in jolts | The rotating blades are not sharp enough or are too rounded on the outside toward the wall | Sharpen or replace the blades, especially if very rounded |
| | Excessive load of waste | Complete the cycle. In the future, stick to a lighter load |
| The time required to cool the treated material is excessively long | The nozzles that send water to the cell are blocked | Clean the nozzles |

| PROBLEM | PROBABLE CAUSES | SOLUTION |
|---|---|---|
| | EV0 and EV1 valves damaged | Check valves functionality (MAINTENANCE SPECIALIST) |
| | Inconsistency or lack of water coming from the water supply | Check that there is pressure in the water supply and that there are no shutoff valves closed |
| The treated material has a bad smell such as organic type | Both carbon filters and absolute filter need to be changed | Replace both filters inside the filters group (MAINTENANCE SPECIALIST). |
| The treated material has a strong smell of toast | The material was kept at too high temperature (approx. 200 ° C) due to dysfunction in the system to measure the temperature | Check the status of the temperature sensor |
| Loss of vapors from the cell | Clogging of the duct of the vapors that connects the cell to the absorption columns | Clean the duct. |
| | Filter group fan not working | Restore the fan functionality (MAINTENANCE SPECIALIST) |
| | Both carbon filters and absolute filter are clogged | Replace both filters inside the filters group (MAINTENANCE SPECIALIST). |
| | The lid gasket is damaged. | Replace the gasket |
| Loss of vapors from the cell and loss of water under the machine | The water level in the first absorption column has exceeded the overflow level and reached the height of the duct of the vapors; in these situations, there is no suction | Check that there are no blockages in the column. Remove sediment on the bottom of the tank; decrease, if excessive, the flow of water sent to the column (MAINTENANCE SPECIALIST) |
| | The circular filter is clogged | Clean the circular filter |
| | The conduit of the vapors that connects the cell to the absorption column is clogged | Open the lid and remove the deposits in the vapors conduct |
| | The filter fan is not working properly or not working at all, clogged filters | Check the filter fan and the connection with the absorption column; check the air filter and replace it if necessary (MAINTENANCE SPECIALIST). |
| It seems that the machine heats less than usual | To much water in the waste load or waste too light | Prepare the load more evenly |
| | Loss of water in the cell for failure of the EV1 and EV0 valves. | Remove and clean the valves (MAINTENANCE SPECIALIST). |
| | Both carbon filters and absolute filter are clogged | Replace both filters inside the filters group |

| PROBLEM | PROBABLE CAUSES | SOLUTION |
|----------------|--|---|
| | The electric heaters that heat the walls of the loading mouth do not work. | Check the fuses and replace the faulty resistors. |