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newster. **NW5**

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HOSPITAL
SOLID
WASTE

 newster.[®]
The new age sterilization

www.newstergroup.com

newster. **NW5**Guarda il video 

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The sterilizing unit is designed for on-site installation up to 150 beds hospitals.

The sturdy design allows the sterilizer to work for several shifts, each lasting 35-40 minutes. The power inverter allows to reduce energy consumption. The automatic treatment takes place in the vessel, which is hermetically sealed during the cycle. The specially designed discharging system allows the operator to easily manage the treated waste without risks.


The unit can be equipped with an optional water recycling system, which increases effectiveness besides further reducing running costs.

Up to 150
beds hospital

newster® NW5 NW15 NW50

Cycle process

- 1. Air inlet
- 2. Temperature sensor
- 3. Water inlet
- 4. Sterilization cell
- 5. Grinder
- 6. Vapour
- 7. Water outlet
- 8. Vapour
- 9. Air
- 10. Air outlet
- 11. Absolute filter
- 12. Carbon activated filter

Technical Specifications

newster. NW5

Sterilization method	Frictional Heat Treatment (Unpressurized moist heat)
Heating method	by impact and friction of the waste
Processing potential	15 kg/h ¹ 240 lt/h
External aspect of treated waste	Homogenous small-sized granules
Final volume of treated waste	20-25% of initial volume
Final weight of treated waste	70-75% of initial weight
Sterilization vessel volume	Roughly 100 liters Ø = 480 mm, H = 550 mm
Overall control system	Programmable logic controller (PLC)
Temperature measurement	By variable-resistance sensors
Cycle recording	Time-temperature flow recording
Final cooling	Roughly 90C° through H ₂ O evaporation
Dust abatement	In humid environment
Vessel lid closure	Manual, mechanical blocking and double safety device
Safety devices	Mechanical blocking system of sterilization vessel lid; power supply to main engine is cut off in case of lid blocking failure. Low voltage command and control panel with automatic power cut-off if panels are open; electrical resistor heating system to be used in case of emergency stop
Rotating blades and fixed contrblades	Made of special metal alloys
H₂O Consumption	Roughly 60 lt/h Roughly 50 lt/day with water recycling system (optional)
Water discharge (diameter sewer pipe equipped with a trap)	40 mm
Power consumption	overall max 20 kW average 13 kW/h

Size and weight	Sterilizer	80 x 160 x 130 cm
	Filter group	80 x 35 x 110 cm
	Electrical board	60 x 45 x 120 cm
	Control panel	-
	Master switch	-
	Total weight	740 kg

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The sterilizing unit is designed for on-site installation up to 300 beds hospitals.

The sturdy design allows the sterilizer to work for several shifts, each lasting 35-40 minutes. The power inverter allows to reduce energy consumption. The automatic treatment takes place in the vessel, which is hermetically sealed during the cycle.


The specially designed discharging system allows the operator to easily manage the treated waste without risks. The unit can be equipped with an optional water recycling system, which increases effectiveness besides further reducing running costs.

Up to 300
beds hospital

newster® NW5 NW15 NW50

Cycle process

- 1. Air inlet
- 2. Temperature sensor
- 3. Water inlet
- 4. Sterilization cell
- 5. Grinder
- 6. Vapour
- 7. Water outlet
- 8. Vapour
- 9. Air
- 10. Air outlet
- 11. Absolute filter
- 12. Carbon activated filter

Technical Specifications

newster. NW15

Sterilization method	Frictional Heat Treatment (Unpressurized moist heat)												
Heating method	by impact and friction of the waste												
Processing potential	30/40 kg/h ¹ 340 lt/h												
External aspect of treated waste	Homogenous small-sized granules												
Final volume of treated waste	20-25% of initial volume												
Final weight of treated waste	70-75% of initial weight												
Sterilization vessel volume	Roughly 170 liters Ø = 500 mm, H = 856 mm												
Overall control system	Programmable logic controller (PLC)												
Temperature measurement	By variable-resistance sensors												
Cycle recording	Time-temperature flow recording												
Final cooling	Roughly 95C° through H ₂ O evaporation												
Dust abatement	In humid environment												
Vessel lid closure	Manual, mechanical blocking and double safety device												
Safety devices	Mechanical blocking system of sterilization vessel lid; power supply to main engine is cut off in case of lid blocking failure. Low voltage command and control panel with automatic power cut-off if panels are open; electrical resistor heating system to be used in case of emergency stop												
Rotating blades and fixed contrblades	Made of special metal alloys												
H₂O Consumption	Roughly 75 lt/h Roughly 50 lt/day with water recycling system (optional)												
Water discharge (diameter sewer pipe equipped with a trap)	40 mm												
Power consumption	overall max 30 kW average 18 kW/h												
Size and weight	<table border="0"> <tr> <td>Sterilizer</td> <td>100 x 200 x 160 cm</td> </tr> <tr> <td>Filter group</td> <td>incorporated or separated</td> </tr> <tr> <td>Electrical board</td> <td>50 x 80 x 160 cm</td> </tr> <tr> <td>Control panel</td> <td>-</td> </tr> <tr> <td>Master switch</td> <td>-</td> </tr> <tr> <td>Total weight</td> <td>1060 kg</td> </tr> </table>	Sterilizer	100 x 200 x 160 cm	Filter group	incorporated or separated	Electrical board	50 x 80 x 160 cm	Control panel	-	Master switch	-	Total weight	1060 kg
Sterilizer	100 x 200 x 160 cm												
Filter group	incorporated or separated												
Electrical board	50 x 80 x 160 cm												
Control panel	-												
Master switch	-												
Total weight	1060 kg												

¹ - Depending on the percentage of humidity and density

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The sterilizing unit is designed for on-site installation up to 400 beds hospitals and external centralized plants.

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
The unit can be equipped with an optional water recycling system, which increases effectiveness besides further reducing running costs. **The system can also be equipped with a platform lift.**

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
Technical Specifications

newster. NW30

Sterilization method	Frictional Heat Treatment (Unpressurized moist heat)												
Heating method	by impact and friction of the waste												
Processing potential	70-80 kg/h ¹ 700 lt/h												
External aspect of treated waste	Homogenous small-sized granules												
Final volume of treated waste	20-25% of initial volume												
Final weight of treated waste	70-75% of initial weight												
Sterilization vessel volume	Roughly 400 liters Ø = 740 mm, H = 900 mm												
Overall control system	Programmable logic controller (PLC)												
Temperature measurement	By variable-resistance sensors												
Cycle recording	Time-temperature flow recording												
Final cooling	Roughly 95C° through H ₂ O evaporation												
Dust abatement	In humid environment												
Vessel lid closure	Manual, mechanical blocking and double safety device												
Safety devices	Mechanical blocking system of sterilization vessel lid; power supply to main engine is cut off in case of lid blocking failure. Low voltage command and control panel with automatic power cut-off if panels are open; electrical resistor heating system to be used in case of emergency stop												
Rotating blades and fixed contrblades	Made of special metal alloys												
H2O Consumption	Roughly 50 lt/day with water recycling system												
Water discharge (diameter sewer pipe equipped with a trap)	50 mm												
Power consumption	overall max 80 kW average 45-50 kW/h												
Size	<table border="0"> <tr> <td>Sterilizer</td> <td>140 x 250 x 180 cm</td> </tr> <tr> <td>Filter group</td> <td>80 x 50 x 110 cm</td> </tr> <tr> <td>Electrical board</td> <td>-</td> </tr> <tr> <td>Control panel</td> <td>80 x 45 x 110 cm</td> </tr> <tr> <td>Master switch</td> <td>50 x 80 x 180 cm</td> </tr> <tr> <td>Total weight</td> <td>2200 kg</td> </tr> </table>	Sterilizer	140 x 250 x 180 cm	Filter group	80 x 50 x 110 cm	Electrical board	-	Control panel	80 x 45 x 110 cm	Master switch	50 x 80 x 180 cm	Total weight	2200 kg
Sterilizer	140 x 250 x 180 cm												
Filter group	80 x 50 x 110 cm												
Electrical board	-												
Control panel	80 x 45 x 110 cm												
Master switch	50 x 80 x 180 cm												
Total weight	2200 kg												

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DATA SHEET

newster® NW50

Sterilizer for Hospital
Solid Waste

HOSPITAL
SOLID
WASTE

Play video 

newster® **NW50**

Up to 600
beds hospital
and external
centralized
plants

Sterilizer for Hospital Solid Waste


The sterilizing unit is designed for on-site installation up to 600 beds hospitals and external centralized plants.

The sturdy design allows the sterilizer to work for several shifts, each lasting 35-40 minutes. The power inverter allows to reduce energy consumption. The specially designed discharging system allows the operator to easily manage the treated waste without risks.

The unit can be equipped with an optional water recycling system, which increases effectiveness besides further reducing running costs. The system can also be equipped with a **platform lift (see Accessories pag. 13)**.


For larger healthcare facilities Newster Twin, which consists of two mirrored NW50 units sharing a single platform lift, can be supplied.

newster® NW5 NW15 NW50

Cycle process

- 1. Air inlet
- 2. Temperature sensor
- 3. Water inlet
- 4. Sterilization cell
- 5. Grinder
- 6. Vapour
- 7. Water outlet
- 8. Vapour
- 9. Air
- 10. Air outlet
- 11. Absolute filter
- 12. Carbon activated filter

Technical Specifications

newster. NW50

Sterilization method	Frictional Heat Treatment (Unpressurized moist heat)												
Heating method	by impact and friction of the waste												
Processing potential	90-110 kg/h ¹ * 840 lt/h												
External aspect of treated waste	Homogenous small-sized granules												
Final volume of treated waste	20-25% of initial volume												
Final weight of treated waste	70-75% of initial weight												
Sterilization vessel volume	Roughly 460 liters Ø = 806 mm, H = 965 mm												
Overall control system	Programmable logic controller (PLC)												
Temperature measurement	By variable-resistance sensors												
Cycle recording	Time-temperature flow recording												
Final cooling	Roughly 95C° through H ₂ O evaporation												
Dust abatement	In humid environment												
Vessel lid closure	Manual, mechanical blocking and double safety device												
Safety devices	Mechanical blocking system of sterilization vessel lid; power supply to main engine is cut off in case of lid blocking failure. Low voltage command and control panel with automatic power cut-off if panels are open; electrical resistor heating system to be used in case of emergency stop												
Rotating blades and fixed contrblades	Made of special metal alloys												
H₂O Consumption	Roughly 50 lt/day with waterrecycling system												
Water discharge (diameter sewer pipe equipped with a trap)	50 mm												
Power consumption	overall max 90 kW average 55-65 kW/h												
Size and weight	<table border="0"> <tr> <td>Sterilizer</td> <td>120 x 250 x 280 cm</td> </tr> <tr> <td>Filter group</td> <td>80 x 35 x 160 cm</td> </tr> <tr> <td>Electrical board</td> <td>-</td> </tr> <tr> <td>Control panel</td> <td>80 x 45 x 110 cm</td> </tr> <tr> <td>Master switch</td> <td>50 x 80 x 180 cm</td> </tr> <tr> <td>Total weight</td> <td>2700 kg</td> </tr> </table>	Sterilizer	120 x 250 x 280 cm	Filter group	80 x 35 x 160 cm	Electrical board	-	Control panel	80 x 45 x 110 cm	Master switch	50 x 80 x 180 cm	Total weight	2700 kg
Sterilizer	120 x 250 x 280 cm												
Filter group	80 x 35 x 160 cm												
Electrical board	-												
Control panel	80 x 45 x 110 cm												
Master switch	50 x 80 x 180 cm												
Total weight	2700 kg												

1 - Depending on the percentage of humidity and density

*Excluding handling

The above specifications are intended for information purposes only and may be subject to change without prior notice.

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

DATA SHEET

newster. **NW50S**

Sterilizer for Hospital
solid waste

HOSPITAL
SOLID
WASTE

 **newster.**
The new age sterilization

www.newstergroup.com

newster. **NW50S**

Sterilizer for Hospital solid waste

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

The unit can be equipped with an optional water recycling system, which increases effectiveness besides further reducing running costs. **The system can also be equipped with a platform lift (see Accessories pag. 13).** For larger healthcare facilities Newster Twin, which consists of two mirrored NW50 units sharing a single platform lift, can be supplied.

**Up to 600
beds hospital
and external
centralized
plants**

newster. NW5 NW15 NW50 Cycle process

- | | | |
|-----------------------|-----------------|-----------------------------|
| 1. Air inlet | 5. Grinder | 9. Air |
| 2. Temperature sensor | 6. Vapour | 10. Air outlet |
| 3. Water inlet | 7. Water outlet | 11. Absolute filter |
| 4. Sterilization cell | 8. Vapour | 12. Carbon activated filter |

Technical Specifications

newster. NW50S

Sterilization method	Frictional Heat Treatment (Unpressurized moist heat)												
Heating method	by impact and friction of the waste												
Processing potential	90-110 kg/h ^{1*} 840 lt/h												
External aspect of treated waste	Homogenous small-sized granules												
Final volume of treated waste	20-25% of initial volume												
Final weight of treated waste	70-75% of initial weight												
Sterilization vessel volume	Roughly 460 liters Ø = 806 mm, H = 965 mm												
Overall control system	Programmable logic controller (PLC)												
Temperature measurement	By variable-resistance sensors												
Cycle recording	Time-temperature flow recording												
Final cooling	Roughly 95C° through H ₂ O evaporation												
Dust abatement	In humid environment												
Vessel lid closure	Manual, mechanical blocking and double safety device												
Safety devices	Mechanical blocking system of sterilization vessel lid; power supply to main engine is cut off in case of lid blocking failure. Low voltage command and control panel with automatic power cut-off if panels are open; electrical resistor heating system to be used in case of emergency stop												
Rotating blades and fixed contrblades	Made of special metal alloys												
H2O Consumption	Roughly 50 lt/day with water recycling system												
Water discharge (diameter sewer pipe equipped with a trap)	50 mm												
Power consumption	overall max 90 kW average 55-65 kW/h												
Size	<table border="0"> <tr> <td>Sterilizer</td> <td>120 x 250 x 280 cm</td> </tr> <tr> <td>Filter group</td> <td>80 x 35 x 160 cm</td> </tr> <tr> <td>Electrical board</td> <td>-</td> </tr> <tr> <td>Control panel</td> <td>80 x 45 x 110 cm</td> </tr> <tr> <td>Master switch</td> <td>50 x 80 x 180 cm</td> </tr> <tr> <td>Total weight</td> <td>2700 kg</td> </tr> </table>	Sterilizer	120 x 250 x 280 cm	Filter group	80 x 35 x 160 cm	Electrical board	-	Control panel	80 x 45 x 110 cm	Master switch	50 x 80 x 180 cm	Total weight	2700 kg
Sterilizer	120 x 250 x 280 cm												
Filter group	80 x 35 x 160 cm												
Electrical board	-												
Control panel	80 x 45 x 110 cm												
Master switch	50 x 80 x 180 cm												
Total weight	2700 kg												

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Newster System S.r.l.